

# CAPP

*the Ultimate in pipetting comfort.....!*

# CAPP

## CappAero™ Multi



- *Single channel accuracy and precision*
- *Superior tip mounting and ejection*
- *Robust and fully autoclavable pipette body*

### Difficulties mounting or ejecting tips?

Many multi channel pipette users have problems with tip pick-up and badly mounted tips are the main cause of inaccuracy and channel-to-channel variability. Often, users must hammer the pipette into the tips to force the tips on - which may increase the risk of RSI (repetitive strain injury).



Working with a CappAero Multi, you only have to press the pipette gently into the tips as the double O-rings on the tip bases make the pipette pick up every tip in the row of the rack. You get high accuracy and even volumes across all channels every time - and ejection is a matter of a light push on the levered ejector arm.

### Ensure you get reliable results

The limits for inaccuracy and imprecision for multi channel pipettes normally are double the limits for single channels\*. With CappAero Multi, you get single channel accuracy and precision across all 8,12 or 16 channels due to the special and sturdy construction of the piston assembly.

Often, tips of varying quality are used for economical reasons. As we mentioned, a bad fit between tip and tip base cause inaccuracies.



The double O-rings of the CappAero Multi increase your chances of getting a good seal, even with a tip that is slightly out of shape. However, we do recommend that you use tips of a high consistent quality such as the original CappExpel™ tips that the pipettes are certified with\*.

### Pipettes that are built to last

Ever dropped a multi-channel pipette and damaged a tip base? For most multis, often this means a costly repair and lack of the pipette for an extended period. Tip bases on the Capp Multi pipettes are made of metal and so they are not prone to break or bend.

Should it happen, an individual tip can very easily be exchanged at your bench. The cost is low, and the down time is negligible.

CappAero multi Pipettes are precision instruments made with a great accuracy, and with critical tolerances down to 1/100 of a millimetre. By using very strong materials, and with great care and workmanship, we build the Capp multi channel pipettes to endure the rough conditions of a working day in the laboratory - when they are new, but also day after day, year after year.

You will have much fewer frustrations due to technical problems and also get the lowest lifetime costs of any multi channel pipette in the world...all backed up with a three year warranty.

\*Please see ISO8655-2



3 Year Warranty