

Pipettors

CappAero™ Multi Pipettes

8-channel pipettes				Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no	Standard	Transparent	SoftLine			
13531/01	13531/05	13531/09	Pipette	White ●	2.50/1.16/1.00	1.80/0.64/0.50
13531/02	13531/06	13531/10	8-channel pipette 0.5-10µl	Red ●	1.50/0.80/0.70	0.80/0.35/0.30
13531/03	13531/07	13531/11	8-channel pipette 0.5-50µl	Yellow ●	0.80/0.72/0.60	0.30/0.22/0.20
13531/04	13531/08	13531/12	8-channel pipette 20-200µl	Black ●	0.70/0.70/0.70	0.30/0.18/0.15
13531/04	13531/08	13531/12	8-channel pipette 50-300µl	Black ●	0.70/0.70/0.70	0.30/0.18/0.15
12-channel pipettes				Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no	Standard	Transparent	SoftLine			
13531/20	13531/24	13531/28	Pipette	White ●	2.50/1.16/1.00	1.80/0.64/0.50
13531/21	13531/25	13531/29	12-channel pipette 0.5-10µl	Red ●	1.50/0.80/0.70	0.80/0.35/0.30
13531/22	13531/26	13531/30	12-channel pipette 0.5-50µl	Yellow ●	0.80/0.72/0.60	0.30/0.22/0.20
13531/23	13531/27	13531/31	12-channel pipette 20-200µl	Black ●	0.70/0.70/0.70	0.30/0.18/0.15
13531/23	13531/27	13531/31	12-channel pipette 50-300µl	Black ●	0.70/0.70/0.70	0.30/0.18/0.15
16-channel pipettes				Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no	Standard	Transparent	SoftLine			
13531/40	13531/42	13531/44	Pipette	Violet ●	12.00/3.10/2.00	6.00/1.60/1.00
13531/41	13531/43	13531/45	16-channel pipette 0.2 - 2µl	Red ●	1.50/0.80/0.70	0.80/0.35/0.30
13531/41	13531/43	13531/45	16-channel pipette 5 - 50µl	Red ●	1.50/0.80/0.70	0.80/0.35/0.30
6-channel HLA pipettes				Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no	Standard	Transparent	SoftLine			
13531/50	13531/51	13531/52	Pipette	White ●	2.50/1.16/1.00	1.80/0.64/0.50
13531/50	13531/51	13531/52	6-channel pipette 0.5 - 10µl	White ●	2.50/1.16/1.00	1.80/0.64/0.50

*Maximum random error (SISO8655) and **Maximum systematic error (ISO8655) at minimum, middle and maximum volume

CappExpel™ Pipette tips for 16 channel pipettes

Product no	Description	Volumes	Size	Packaging
13531/53	CappExpel 02/384 - 0.2 to 30µl (length 45 mm)	● ●	1 x 384	Sterile, racks w/loose lid
13531/54	CappExpel 02/384 - 0.2 to 30µl (length 45 mm)	● ●	10 x 384	Sterile, racks w/loose lid
13531/55	CappExpel 10/384 - 0.5 to 60µl (length 50 mm)	● ●	1 x 384	Sterile, racks w/loose lid
13531/56	CappExpel 10/384 - 0.5 to 60µl (length 50 mm)	● ●	10 x 384	Sterile, racks w/loose lid

CappWash™ - 16 channel pipettes

Product no	Description
13531/57	16-channel washer for 384-well microplates

Yorlab pipette, calibration and service

Please call for a quote