

Eyeshields dispenser

Ideal storage and dispensing box for visitors eyeshields, keeping them readily available and dust free. These compact dispensers are made from smoked acrylic in two sizes. The small dispenser holds 5 pairs and measures 190 x 170 x 90mm (HxWxD). Alternatively, the larger dispenser holds 10 pairs and measures 320 x 170 x 90mm, respectively. Either size is suitable for use with "Coverspecs" and other eyeshields.

Product No	Size
10816/01	Large
10816/02	Small

Face Shields

Extremely light, fully adjustable head harness and browguard made from a mix of fibreglass, nylon and PVC, giving durability and flexibility. Visors are all 230mm (9") deep. Choose from polycarbonate for impact resistance (to BS EN1662B39) or acetate for protection from most chemicals (to BS EN1662F3)



Product No	Description
10813/01	230mm clear polycarbonate faceshield
10813/03	230mm clear acetate faceshield
10813/02	Spare polycarbonate visor
10813/11	Spare acetate visor

Heat Resistant Gloves

Textured Palm

Product No 10817/01
Made from cut-resistant Kevlar, these gloves will withstand temperatures up to 100°C. A textured palm gives a superior grip.



Textured Palm & Cotton Liner

Product No 10818/01
Polysafe gloves with similar qualities to the above, but a cotton liner for wearer comfort and greater protection from heat, up to 250°C.



Heat & Water Resistant

Product No 10819/01
Polysafe Kevlar gloves with similar qualities to 10818/01 with cut and heat resistance to 250°C and with nitrile interlining providing resistance to water and many chemicals. A cotton liner improves wearer comfort.

Product No 10820/01

Heat resistant gloves offering protection up to 900°C.



Anti Static Gloves

Lint-free Nylon

Washable gloves, made of lightweight 100% stretch nylon offering both comfort and protection

Product No	Size
11370/01	7 1/2 (Ladies)
11370/02	8 1/2 (Men)

Cryogenic Gloves

Cotton lined, polyolefin/polyester micro-fibre insulated gloves. These gloves are highly flexible even at temperatures to -190°C, are washable, and water resistant. Supplied by the pair in the sizes indicated below.

Product No	Size
10490/01	Small
10490/02	Medium
10490/03	Large

