

Safety

High Neck Labcoats

Unisex

White coats in polycotton with off-set front stud fastening and stretch cuffs. Includes two hip pockets and one breast pocket. 110cm long

Product No	Size-cm
11342/01	80
11342/02	84
11342/03	88
11342/04	92
11342/05	96
11342/06	100
11342/07	104
11342/08	108
11342/09	112
11342/10	116
11342/11	124
11342/12	132
11342/13	140

Headwear

Unisex

A choice of white hats suitable for use in food preparation areas for safety and hygiene.



Product No	Size
11343/01	Coverall with net top and polyester/cotton peak
11343/02	Polyester/cotton cap with snood

Wellington Boots

Non-safety

White, high grade PVC nitrile boots resistant to petrol, oils, fats, chemicals, acids, etc. A smooth finish makes them easy to clean.



Mens		Ladies	
Product No	Size	Product No	Size
11344/01	6	11344/08	3
11344/02	7	11344/09	4
11344/03	8	11344/10	5
11344/04	9	11344/11	6
11344/05	10	11344/12	7
11344/06	11		
11344/07	12		

Overshoes

Elasticated, hardwearing PVC overshoes, suitable for clean or wet areas. Supplied in packs of 100



Product No	Size-inches
10829/01	14
10829/02	16

Sleeve Protectors

Product No 11345/01

Disposable blue polythene sleeves to protect forearms from occasional splashes. Supplied in packs of 2000



Dust/Mist Respirator

Product No 11347/01

An effective, CE marked, lightweight, respirator which provides protection against fine respirable dusts and water-based mists up to 10 times the Occupational Exposure Limit (APF). Its convex shape, twin fitting straps, foam nose seal with soft aluminium nose clip, ensure comfort and a good seal. Conforms to EN149FFP2S. Supplied in packs of 20 masks

