

Sieving

Measuring Scoops

Rigid Polypropylene scoops.

| Product No | Capacity-ml |
|------------|-------------|
| 10831/01 | 10 |
| 10831/02 | 25 |
| 10831/03 | 50 |
| 10831/04 | 100 |
| 10831/05 | 250 |
| 10831/06 | 500 |
| 10831/07 | 1000 |

Test Sieves

Stainless steel framed, 200mm diameter, nesting sieves with precision woven wire mesh and conforming to BS410:86. Supplied individually boxed with certificate of conformity. We can also supply sieves with stainless steel frames and other mesh sizes



| Product No | Mesh size-mm | Product No | Mesh size-µm |
|------------|--------------|------------|--------------|
| 10861/01 | 16.00 | 10862/01 | 850 |
| 10861/02 | 13.20 | 10862/02 | 710 |
| 10861/03 | 11.20 | 10862/03 | 600 |
| 10861/04 | 9.50 | 10862/04 | 500 |
| 10861/05 | 8.00 | 10862/05 | 425 |
| 10861/06 | 6.70 | 10862/06 | 355 |
| 10861/07 | 5.60 | 10862/07 | 300 |
| 10861/08 | 4.75 | 10862/08 | 250 |
| 10861/09 | 4.00 | 10862/09 | 212 |
| 10861/10 | 3.35 | 10862/10 | 180 |
| 10861/11 | 2.80 | 10862/11 | 150 |
| 10861/12 | 2.36 | 10862/12 | 125 |
| 10861/13 | 2.00 | 10862/13 | 106 |
| 10861/14 | 1.70 | 10862/14 | 90 |
| 10861/15 | 1.40 | 10862/15 | 75 |
| 10861/16 | 1.18 | 10862/16 | 63 |
| 10861/17 | 1.00 | 10862/17 | 53 |
| | | 10862/18 | 45 |
| | | 10862/19 | 38 |
| | | 10862/20 | 32 |
| | | 10862/21 | 25 |
| | | 10862/22 | 20 |

Lids & Receivers

Brass lids and receivers for 200mm sieves.

| Product No | Description |
|------------|-----------------------|
| 10863/01 | Lid |
| 10863/02 | Bottom receiver |
| 10863/03 | Intermediate receiver |

Sieve Shakers

Suitable for any sieves up to 315mm diameter, an integral timer can be used to ensure sieving continues for a set period. The benchtop shaker will hold ten 200mm diameter sieves or six 300mm sieves, in addition to a lid and receiver. It has a random, multi-directional motion to sieve granules up to 20mm. The floor-standing model will hold eleven 200mm diameter sieves or seven 300mm sieves. The sieving action is both vertical and horizontal.



| Product No | Description |
|------------|----------------|
| 10969/01 | Benchtop |
| 10969/02 | Floor standing |

Sieve Brush



| Product No | Description |
|------------|--|
| 10129/01 | Double ended sieve brush with soft and stiff bristles for use with both fine and coarse sieves respectively. |