

Spectrophotometer Accessories

Sample Holders

Product No	For Series	Description
11375/11	All	Single 100mm cell holder
11375/12	All	Rectangular cell holder (10 to 100mm cells)
11375/13	All	Cylindrical cell holder (10 to 100mm)
11375/14	All	Holder for solid sample or films
11375/15	All	Ultra-microcell holder: 50µl or less (Requires cell 11375/09)
11375/16	All	SIP system flowcell holder for 10mm cells
11375/17	1000	Holder for up to 50mm flowcells
11375/18	1000	Ultra microcell holder for samples up to 50µl

Cell Changers & Holders

Product No	For Series	Description
11375/21	All	Manual, 4 microcells (with 4 x 10mm cell holder)
11375/22	All	Manual, 3 x 10-40mm cells with holder for water analysis
11375/23	1000	Manual, 4 cells (requires holder)
11375/24	1000	Holder (4 x 10mm cells) for 11375/23
11375/25	1000	Holder (up to 4 x 50mm cells) for 11375/23
11375/26	2000, 3000	Automatic, 4 cells (requires holder)
11375/27	2000, 3000	Holder (4 x 10mm cells) for 11375/26
11375/28	2000, 3000	Holder (4 x 50mm cells) for 11375/26
11375/29	2000, 3000	Automatic, 6 cells (with 10mm cell holder)
11375/30	2000, 3000	Automatic, 8 x 10mm cells with lid & manifold
11375/31	2000, 3000	Sample compartment lid and pipe guide for tablet dissolution (4 & 6 cell changers)

Temperature Control

With water thermostating, accessory 11375/66 will also be required.

Product No	For Series	Description
11375/61	All	Holder for single, 10mm SIP cell
11375/62	All	Holder for manual 4 x 10mm cell changer
11375/63	2000, 3000	Holder for automatic 4 x 10mm cell changer
11375/64	2000, 3000	Peltier controller and 10mm cell holder
11375/65	2000, 3000	Peltier controller & auto change
11375/66	All	Cell compartment lid and water pipe port

Spectrophotometer Accessories

SIP & Sipette and Batch Sampling

Product No	For Series	Description
11375/41	2000	Batch sampler for 40 samples, with pump
11375/42	3000	Batch sampler for 40 samples, with pump
11375/43	2000	Batch sampler for 100 samples, with pump
11375/44	2000, 3000	Batch sampler for 100 samples, with pump
11375/45	2000, 3000	Return Sipette system (requires cell & holder)
11375/46	2000	As above, with Viton components
11375/47	1000, 2000	Holder for 10mm SIP/Sipette cell
11375/48	2000, 3000	Holder for up to 50mm flowcells
11375/49	All	Holder for SIP/Sipette cell up to 50mm
11375/50	1000	Tube kit for SIP cell system (2 sets)
11375/51	2011, 2031	Tube kit for Sipette cell system (2 sets)
11375/52	2021, 2041, 3000	Tube kit for Sipette cell system (2 sets)
11375/53	2000, 3000	Viton tube kit (2 sets)
11375/54	2000, 3000	10mm silica Micro Flow Through Cell
11375/55	2000, 3000	10mm silica Sipette cell, piping and connectors
11375/56	2000, 3000	20mm silica Sipette cell, piping and connectors
11375/57	2000, 3000	40mm silica Sipette cell, piping and connectors

Sample Stirring

Product No	For Series	Description
11375/71	All	Electronic controller for sample stirring
11375/72	All	Cell holder for sample stirring
11375/73	All	Water thermostatted holder for stirring

Lamps

Product No	For Series	Description
11375/81	All	Deuterium lamp with hours-used indicator
11375/82	All	Tungsten Halide lamp (pack of 2)