


Stands, Scaffolds and Racks

Laboratory Jacks

Hand operated jacks of chemical resistant, stainless steel.



Product No	Plate area-mm	Min/max height- mm	Max load-kg
10839/01	100 x 120	55 / 210	50
10839/02	130 x 160	55 / 270	70
10839/03	180 x 220	73 / 375	90

Retort Stands

Rectangular Base

The base plate measures 200 x 315mm and the rod is 700mm tall, 16mm diameter. Total weight is 5.2kg. Suitable for light duty.



H-Frame Base

The foot measures 430 x 415mm and the rod is 700mm tall, diameter 25mm. Total weight is 5.8kg. Suitable for medium/heavy duty.



Product No	Description
10711/01	Rectangular base & rod
10711/02	H-Frame base & rod
10711/03	Clamp for rods 13 to 25mm dia
10711/04	D cross sleeve for rods 13 to 16mm dia

Retort Stands & Supports

Made of steel with rubber feet and finished in blue acrylic. A 10mm threaded hole is located midway along the short side



Product No	Dimensions	Weight
10244/01	115 x 75mm	0.30kg
10244/02	160 x 100mm	0.75kg
10244/03	200 x 125mm	1.25kg
10244/04	250 x 160mm	2.00kg
10244/05	315 x 200mm	3.00kg

Tripod Bases

Manufactured in cast iron with rubber feet and finished in blue acrylic. A 10mm threaded hole is located in the centre



Rectangular Bases

Product No	Dimensions	Weight
11358/01	115 x 75mm	0.30kg
11358/02	160 x 100mm	0.75kg
11358/03	200 x 125mm	1.25kg
11358/04	250 x 160mm	2.00kg
11358/05	315 x 200mm	3.00kg

Rods

With 10mm screwthread, suitable for bases on earlier page

Rods

Product No	Material	Diameter x length-mm
10246/01	Mild steel	10 x 500
10246/02	Mild steel	12.5 x 600
10246/03	Mild steel	12.05 x 760
10246/04	Mild steel	12.5 x 1000
10246/05	Stainless steel	10 x 500
10246/06	Stainless steel	12.5 x 600
10246/07	Stainless steel	12.05 x 760
10246/08	Stainless steel	12.5 x 1000