

# Stands, Scaffolds and Racks

## Retort Rods

Zinc plated steel rods, rounded both ends, which may be linked together with bossheads (to be ordered separately) to form a secure scaffold frame.

Product No	Diameter x length-mm
10247/01	10 x 500
10247/02	12.5 x 600
10247/03	12.5 x 750
10247/04	12.5 x 1000

## Bosshead, Straight



**Product No** 10241/01  
Made of natural finished alloy, with plastic headed hand screws. Will hold rods up to 16mm diameter

## Retort Clamps

### General Purpose, Four-Prong

**Product No** 10242/01  
Made of natural alloy with cork-lined prongs, having a maximum jaw opening of 80mm. Rod diameter is 8mm.



### General Purpose, Three-Prong

**Product No** 10243/01  
Made of natural alloy with rubber-covered prongs, having a maximum jaw opening of 70mm. Rod diameter is 8mm



## Reaction Vessel Clamp

**Product No** 10249/01  
Made of alloy and zinc-plated steel with heat-resistant covered jaws. Each clamp is supplied with two extension rods. One is for vessels up to 5 litres and the second for up to 20 litres

## Beaker Clamp

**Product No** 11359/01  
Suitable for use with beakers up to 5 litres capacity, the clamp is made of alloy and stainless steel with heat-resistant covered jaws. Will accept items up to 175mm diameter

## Retort Rings, Plain Stem

Mild steel with zinc plated rings. The rods are 140mm long, 8mm diameter

Product No	Ring ID-mm
10250/01	38
10250/02	50
10250/03	56
10250/04	62
10250/05	69
10250/06	75
10250/07	82
10250/08	100
10250/09	106
10250/10	125
10250/11	131
10250/12	150



## Retort Rings, Including Bosshead

As above, with alloy bosshead

Product No	Ring ID-mm
10253/01	50
10253/02	75
10253/03	100
10253/04	125
10253/05	150

