

# Stands, Scaffolds and Racks

## Scaffolding

### Heavy Duty

Made of corrosion resistant stainless steel, the 25mm tubing can provide a robust frame. The optional castors can be fitted to the tubing ends to allow assemblies of glass and other apparatus to be transported with ease and in complete safety



Tubing, 25mm OD

Product No	Length-mm
10731/01	250
10731/02	500
10731/03	750
10731/04	1000
10731/05	1500
10731/06	2000
10731/07	3000
10731/08	4000
10731/09	5000
10731/10	6000

### Vessel Base Supports

Product No	For vessel OD-mm
10731/18	150
10731/19	200
10731/20	250

### Other Fittings

Product No	Description
10731/11	Plastic end cap for tubing
10731/15	Bosshead, unequal (30mm max to 20mm max)
10731/12	Duralim knuckle 25mm bosshead connector
10731/16	Flange support clamp, DN150
10731/17	Flange support clamp, DN200
10731/13	Castor wheel
10731/14	Castor wheel, with brake

## Burette Clamps

A polypropylene assembly that will fit support rods with diameters between 8 - 14mm, and holding one or two burettes, as indicated. Spring loaded grips hold the burette securely while allowing easy removal. Graduations and meniscus remain visible. Suitable retort rods and bases may be found on earlier pages, if required.

Product No	No. of Positions
10835/01	1
10835/02	2

## Funnel Holders

Single or double position polypropylene holders suitable for funnel diameters from 25mm to 170mm. (Will fit retort rods from 8mm to 14mm. Please see earlier page )

Product No	No. of Positions
10833/01	1
10833/02	2



## Separating Funnel Holders

Single position, similar to above, but with open side for separating funnels of 75mm diameter and above. Will fit retort rods from 8 - 14mm, see earlier pages.

Product No	No. of Positions
10834/01	1

