

Stands, Scaffolds and Racks

Test Tube Racks

Nylon Coated Wire


Autoclavable to 130°C. Other sizes, configurations and materials are available



Z-Rack

Product No **10932/01**

Ideal for general use, these aluminium racks accept tubes of various sizes. The 'Z' shape allows nesting for space saving storage. The rack provides five 28mm holes, six 19mm holes and eight 13mm holes.



Product No	Tube Dia-mm	Array	Dimensions (H xW xL)-mm
10345/01	13	4 x 6	50 x 80 x 110
10345/02	13	4 x 20	50 x 80 x 360
10345/03	13	3 x 12	50 x 62 x 220
10345/04	16	3 x 4	100 x 76 x 94
10345/05	16	4 x 6	100 x 98 x 139
10345/06	16	3 x 12	100 x 76 x 272
10345/11	17	5 x 10	85 x 123 x 228
10345/07	20	3 x 6	100 x 88 x 160
10345/08	20	3 x 12	100 x 88 x 310
10345/09	25	3 x 4	100 x 102 x 128

Microtube Racks

Floating

These racks are made of polypropylene and resistant to temperatures from -70°C to +100°C and are autoclavable. The racks will float with a full load of filled tubes. The rack legs raise the tubes above the work surface. Square racks are suited for work in baths, circular racks for work in beakers. Supplied in packs of 4



Traditional Racks

Made of autoclavable polypropylene, these racks are lightweight and virtually unbreakable

Product No	Holes-mm	Height-mm
11399/01	12 x 12	80
11399/02	12 x 20	80
11399/03	24 x 20	95
11399/04	10 x 18	74

Product No	10873/01	10873/02	10873/03	10873/04
Colour	Black	White	White	White
Tube size-ml	0.5	1.0, 1.2, 1.5 & 2.0	1.0, 1.2, 1.5 & 2.0	1.0, 1.2, 1.5 & 2.0
Shape	Square	Square	Circular	Circular
Size-mm	103 x 103 x 65	103 x 103 x 65	Dia 66	Dia 66
Array	4 x 4	4 x 4	8	20

