

Stoppers

Stoppers


Rubber & Silicone Rubber

All stoppers are manufactured to BS 2775. Dimensions are indicated in the following tables:

BS No	Bottom dia-mm	Top dia- mm	Length mm	BS No	Bottom dia-mm	Top dia- mm	Length mm
3	3	4.5	12	33	33	38.5	38.5
4	4	5.8	14	35	35	40.5	38.5
5	5	7	16	37	37	42.5	38.5
6	6	8	16	39	39	44.5	38.5
7	7	9	16	41	41	47	42
8	8	10.5	20	43	43	49	42
9	9	11.5	20	45	45	51	42
10	10	12.5	20	49	49	55	42
11	11	14	24	53	53	59.5	45.5
13	13	16	24	57	57	64	49
15	15	18	24	62	62	69	49
17	17	20.5	26	67	67	74	49
19	19	22.5	28	71	71	78	49
21	21	24.5	28	75	75	85	54
23	23	26.5	28	79	79	89	54
25	25	28.5	28	86	86	96	54
27	27	31	32	92	92	102	54
29	29	33	32	99	99	108	54
31	31	36	35	105	105	117	63.5

Natural Rubber, Solid

Product No	BS No	Pack size
10216/01	3	10
10216/02	4	10
10216/03	5	10
10216/04	6	10
10216/05	7	10
10216/06	8	10
10216/07	9	10
10216/08	10	10
10216/09	11	10
10216/10	13	10
10216/11	15	10
10216/12	17	10
10216/13	19	10
10216/14	21	10
10216/15	23	10
10216/16	25	10
10216/17	27	10
10216/18	29	10
10216/19	31	10
10216/20	33	10
10216/21	35	10
10216/22	37	10
10216/23	39	10




Product No	BS No	Pack size
10216/24	41	5
10216/25	43	5
10216/26	45	5
10216/27	49	5
10216/28	53	5
10216/29	57	5
10216/30	62	5
10216/31	67	1
10216/32	71	1
10216/33	75	1
10216/34	79	1
10216/35	86	1
10216/36	92	1
10216/37	99	1
10216/38	105	1

Stoppers

Natural Rubber With One Hole

Sizes 5 to 8 have a 3mm hole. Other sizes have a 5mm hole

Product No	BS No	Pack size
10217/01	5	10
10217/02	6	10
10217/03	7	10
10217/04	8	10
10217/05	9	10
10217/06	10	10
10217/07	11	10
10217/08	13	10
10217/09	15	10
10217/10	17	10
10217/11	19	10
10217/12	21	10
10217/13	23	10
10217/14	25	10
10217/15	27	10
10217/16	29	10
10217/17	31	10
10217/18	33	10
10217/19	35	10
10217/20	37	10
10217/21	39	10



Product No	BS No	Pack size
10217/22	41	5
10217/23	43	5
10217/24	45	5
10217/25	49	5
10217/26	53	5
10217/27	57	5
10217/28	62	5
10217/29	67	1
10217/30	71	1
10217/31	75	1
10217/32	79	1
10217/33	86	1
10217/34	92	1
10217/35	99	1
10217/36	105	1

Natural Rubber, With Two Holes

Sizes 5 to 8 have two 3mm holes. Other sizes have two 5mm

Product No	BS No	Pack size	Product No	BS No	Pack size
10218/01	15	10	10218/15	43	5
10218/02	17	10	10218/16	45	5
10218/03	19	10	10218/17	49	5
10218/04	21	10	10218/18	53	5
10218/05	23	10	10218/19	57	5
10218/06	25	10	10218/20	62	5
10218/07	27	10	10218/21	67	1
10218/08	29	10	10218/22	71	1
10218/09	31	10	10218/23	75	1
10218/10	33	10	10218/24	79	1
10218/11	35	10	10218/25	86	1
10218/12	37	10	10218/26	92	1
10218/13	39	10	10218/27	99	1
10218/14	41	5	10218/28	105	1