

Stoppers

Silicone Rubber, Solid

Product No BS No Pack size

10219/01	3	10
10219/02	4	10
10219/03	5	10
10219/04	6	10
10219/05	7	10
10219/06	8	10
10219/07	9	10
10219/08	10	10
10219/09	11	10
10219/10	13	10
10219/11	15	10
10219/12	17	10
10219/43	19	5
10219/44	21	5
10219/45	23	5
10219/46	25	2
10219/47	27	2
10219/48	29	2
10219/49	31	2
10219/50	33	2
10219/51	35	2
10219/52	37	2
10219/53	39	2
10219/54	41	2



Product No BS No Pack size

10219/55	43	2
10219/56	45	2
10219/57	49	2
10219/58	53	2
10219/59	57	2
10219/60	62	2
10219/31	67	1
10219/32	71	1
10219/33	75	1
10219/34	79	1
10219/35	86	1
10219/36	92	1
10219/37	99	1
10219/38	105	1

PTFE Stoppers

Machined in pure chemical resistant P.T.F.E to fit standard conical joints. Trapezoidal seals give medium to high vacuum performance to 280°C.

Product No Joint size

10736/01	10/19
10736/02	14/23
10736/03	19/26
10736/04	24/29
10736/05	29/32
10736/06	34/35

Polystop Stoppers

Polypropylene interchangeable stoppers for vessels with conical joints. Supplied in packs of 20.

Product No	Joint size	Product No	Joint size
10377/01	10/13	10377/05	19/17
10377/02	12/14	10377/06	24/20
10377/03	14/15	10377/07	29/32
10377/04	16/16	10377/08	34/35

Nesting Stoppers (GUKO sets)

Made from EDPM rubber, this set of eight gaskets will withstand temperatures up to 150°C.

Product No **10840/01**



Contents all measurements in mm

Gasket No	1	2	3	4	5	6	7	8
Large OD	22	29	36	44	53	63	73	84
Small OD	12	16	22	27	33	43	52	61
Small ID	7	9	15	19	24	33	42	50
Height	17	22	25	30	35	35	40	40

Cork Borers

Nickel-plated brass borers

Product No Set of Hole-mm

11389/01	6	4 - 10
11389/02	12	4 - 18
11389/03	18	4 - 25

Sharpener

Metal, suitable for above borers

Product No

11389/11

