

Tally Counter and Test Papers

Tally Counter

Product No 10954/01

A robust, thumb-operated hand-held counter, with finger ring for a secure grip. Counts from 0-9999 and has a re-set knob



Indicator Papers

pH Type 112

Supplied in splash-proof box.

Product No	Range-pH	Graduations-pH	Pack Quantity
10278/50	0 - 14	0.5 (pH 0.9) 1.0 (pH 9 - 14)	200 strips
10278/51	1 - 11	Universal	Roll (7mm x 5M)

Panaphea Plus - Indicator Sticks

Non-curling plastic strips (6x85mm) with colour change pads.

Product No	Range-pH	Graduations-pH	Pack Quantity
10278/55	0 - 14	1.0	100
10278/56	2 - 9	0.5	100
10278/57	7 - 14	0.5	100

Specialist pH and Test Papers

As indicated below. Supplied in a booklet of 100 strips measuring 10x75mm.

Product No	Description	Determination
10278/60	Litmus Blue	pH 8 to 5
10278/61	Litmus Red	pH 5 to 8
10278/62	Congo Red	pH 3 to 5
10278/63	Phenolphthalein	pH > 8.3
10278/64	Lead Acetate	Hydrogen Sulphide
10278/65	Potassium Iodide (starch)	Nitrides

pH Papers - Type CF

Use for unbuffered or weakly buffered solutions. Boxes of 100 strips each measuring 6x80mm.

Product No	Range-pH	Graduations-pH
10384/01	0 to 14	1.0
10384/02	4.5 to 10	0.5

Type CS

Suitable for coloured solutions. Boxes of 200 strips, each measuring 11x100mm.

Product No	Range-pH	Graduations-pH
10385/01	1.0 - 12.0	1.0
10385/02	1.8 - 3.8	0.2 - 0.3
10385/03	3.8 - 5.5	0.2 - 0.3
10385/04	5.2 - 6.8	0.2 - 0.3
10385/05	6.0 - 8.1	0.3
10385/06	8.09 - 9.7	0.2 - 0.3

Type TC

Three colour indicator. Five metre reel of paper 10mm wide.

Product No	Range-pH	Graduations-pH
10386/01	1.0 - 11.0	1.0

Type FR & NR

Single band, supplied in a convenient pack of 4 dispensers containing 5 metres of paper 6mm wide. Just tear off a suitable length, immerse and compare with the colour chart on the side of the reel.

Product No	Range-pH	Graduations-pH
10278/01	1.0 - 14.0	1.0 - 2.0
10278/02	1.0 - 4.0	0.5 - 1.0
10278/03	4.0 - 6.0	0.5
10278/04	6.0 - 8.0	0.5
10278/05	8.0 - 10.0	0.5

Product No 10387/01

A pack containing two 10278/01 plus one of each of the other 10278 reels.

Specialist Test Papers

Supplied in a pack of four dispenser reels each containing 5 metres of indicator paper 6mm wide.

Product No	Description	Determination
10302/01	Litmus Blue	pH8 to 5
10302/02	Litmus Red	ph5 to 8
10302/03	Congo Red	Acid/alkali
10302/04	Phenolphthalein	pH8.3/10
10302/06	Starch Lodide	Chlorine/oxidising agents

As above, but supplied in a carton of 10 packs of 10 books of 20 leaves:

Product No	Description	Determination
10303/01	Limus Blue	
10303/02	Litmus Red	
10303/04	Lead acetate	
10303/05	Starch Lodide	Chlorine/oxidising agents