

## Cryoscope Tubes

**Product No** 10383/01

2.5ml capacity, glass, two-mark cryoscope tubes manufactured to BS3095. Supplied singly




## Hagberg Tubes

**Product No** 10952/01

Precision bore, weight controlled test tubes of borosilicate glass. Used in the Hagberg test for the determination of moisture in grain by the falling number method

## Test Tubes



### Polystyrene

Single use crystal clear tubes. Supplied in packs of 1000.

## Culture Tubes

Disposable, round bottom rimless tubes of borosilicate glass (EXCEPT 10255/03, soda glass). Supplied in cases of 1000 (each case contains four boxes of 250 tubes per box)

Product No	Length x Dia-mm
10255/01	75 x 10
10255/02	75 x 12
10255/03	100 x 13
10255/04	100 x 16
10255/05	125 x 16
10255/06	150 x 16

## Durham Tubes

Rimless tubes of high quality soda glass with a wall thickness of 0.65mm. Supplied in packs of 250



Product No	Type	Length x Dia-mm
10274/01	Conical	105 x 16
10274/02	Cylindrical	100 x 16
10274/03	Cylindrical	75 x 12
10274/06	Cylindrical	150 x 16
10274/04	Coagulometer	55 x 11.5
10274/05	RIA	70 x 11

### Stoppers

Made of polyethylene and supplied in packs of 1000.

Product No	For Tubes	Colour
11398/01	10274/01, /02	Neutral
11398/02	10274/01, /02	Blue
11398/03	10274/01, /02	Yellow
11398/04	10274/01, /02	Red
11398/05	10274/01, /02	Green
11398/06	10274/03, /04	Neutral
11398/07	10274/03, /04	Blue
11398/08	10274/03, /04	Yellow
11398/09	10274/03, /04	Red
11398/10	10274/03, /04	Green
11398/11	10274/05	Neutral
11398/12	10274/06	White
11398/13	10274/06	Blue
11398/14	10274/06	Yellow
11398/15	10274/06	Black
11398/16	10274/06	Red
11398/17	10274/06	Pink

Product No	Length x Dia-mm
10866/01	25 x 6.5
10866/02	30 x 6.5
10866/03	35 x 7.5
10866/04	50 x 7.5