

Test Tubes

Test Tubes

Neutral Borosilicate Glass


Light wall (0.8mm), rimless tubes suitable for general purpose work

Without Rim:



Product No	Length x Dia	Pack quantity
10867/01	75 x 10mm	100
10867/02	75 x 12mm	100
10867/03	100 x 12mm	100
10867/04	100 x 16mm	100
10867/05	125 X 16mm	100
10867/06	150 x 16mm	100
10867/07	150 x 18mm	100
10867/08	150 x 24mm	50

With Rim:



Product No	Length x Dia	Pack quantity
10868/01	75 x 10mm	100
10868/02	75 x 12mm	100
10868/03	100 x 12mm	100
10868/04	100 x 16mm	100
10868/05	125 x 16mm	100
10868/06	150 x 16mm	100
10868/07	150 x 18mm	100
10868/08	150 x 24mm	50

Test Tubes

Heavy wall for maximum mechanical and thermal shock resistance.

Without Rim:



Product No	Length x OD	Wall thickness	Pack quantity
10325/01	70 x 8mm	0.8-1.0mm	100
10325/02	75 x 10mm	0.8-1.0mm	100
10325/03	100 x 10mm	0.8-1.0mm	100
10325/04	75 x 12mm	0.8-1.0mm	100
10325/05	100 x 12mm	0.8-1.0mm	100
10325/06	130 x 14mm	0.8-1.0mm	100
10325/07	130 x 16mm	1.0-1.2mm	100
10325/08	160 x 16mm	1.0-1.2mm	100
10325/09	180 x 18mm	1.0-1.2mm	100
10325/10	150 x 20mm	1.0-1.2mm	100
10325/11	180 x 20mm	1.0-1.2mm	100
10325/26	150 x 25mm	1.0-1.2mm	80
10325/13	200 x 25mm	1.0-1.2mm	50
10325/14	200 x 30mm	1.0-1.4mm	50

With Rim:

Product No	Length x OD	Wall thickness	Pack quantity
10326/01	70 x 8mm	0.8-1.0mm	100
10326/02	75 x 10mm	0.8-1.0mm	100
10326/03	100 x 10mm	0.8-1.0mm	100
10326/04	75 x 12mm	0.8-1.0mm	100
10326/05	100 x 12mm	0.8-1.0mm	100
10326/06	130 x 14mm	0.8-1.0mm	100
10326/07	130 x 16mm	1.0-1.2mm	100
10326/08	160 x 16mm	1.0-1.2mm	100
10326/09	180 x 18mm	1.0-1.2mm	100
10326/10	150 x 20mm	1.0-1.2mm	100
10326/11	180 x 20mm	1.0-1.2mm	100
10326/12	150 x 25mm	1.0-1.2mm	50
10326/13	200 x 25mm	1.0-1.2mm	50
10326/14	200 x 30mm	1.0-1.4mm	50