

Test Tubes

Extra strong, heavy wall rimless test tubes of soda glass. 10596/02 and 10596/03 are suitable for Wasserman tests. 10596/06 tubes are suitable for Kahn tests

Product No	Length x OD	Wall thickness	Pack quantity
10596/01	50 x 6.5mm	0.65mm	250
10596/02	50 x 10mm	1.0mm	250
10596/03	75 x 10mm	1.0mm	100
10596/04	90 x 10mm	1.0mm	100
10596/05	50 x 12mm	1.0mm	100
10596/06	75 x 12mm	1.0mm	100
10596/07	100 x 12mm	1.0mm	100
10596/08	75 x 16mm	1.0mm	100
10596/09	100 x 16mm	1.0mm	100
10596/10	125 x 16mm	1.0mm	100
10596/11	150 x 16mm	1.0mm	100
10596/12	125 x 18mm	1.0mm	100
10596/13	150 x 18mm	1.0mm	100
10596/14	150 x 24mm	1.2mm	50

Test Tube Brushes

Plain

Nylon bristles with fan shaped head. Available in two sizes indicated below. Supplied in packs of 10.

Product No	For Tube Diameter-mm
10126/01	10 - 12
10126/02	16 - 25

Cotton Tipped

Nylon bristles with cotton tip on wire handle available in two sizes. Supplied in packs of 10.

Product No	For Tube Diameter-mm
10127/01	12 - 16
10127/02	19 - 25



Test tube caps

Rust-proof anodised aluminium caps with an internally fitted chrome-nickel spring which securely grips the inside of the tube. Caps provide a sterile (but not hermetic) seal and are colour coded for easy sample identification. Other sizes and colours are available to order. Supplied in packs of 100.



Product No	Colour	Size
10869/01	Silver	12/13mm
10869/02	Silver	15/16mm
10869/03	Silver	16/17mm
10869/04	Silver	17/18mm
10869/05	Silver	19/20mm
10869/06	Red	12/13mm
10869/07	Red	15/16mm
10869/08	Red	16/17mm
10869/09	Red	17/18mm
10869/10	Red	19/20mm
10869/11	Blue	12/13mm
10869/12	Blue	15/16mm
10869/13	Blue	16/17mm
10869/14	Blue	17/18mm
10869/15	Blue	19/20mm
10869/16	Gold	12/13mm
10869/17	Gold	15/16mm
10869/18	Gold	16/17mm
10869/19	Gold	17/18mm
10869/20	Gold	19/20mm

Test Tube Holders

Brass Cross-Over

Simple, durable holders with sprung brass tube grip and smooth wooden handle. Will hold tubes up to approximately 35mm diameter

Product No 10929/01



Wire

Spring type holders of galvanised wire. Will hold tubes up to approximately 50mm diameter.

Product No 10956/01

