


# Thermometers

## Stirring/Immersion Thermometers

### Mercury in Glass

General Purpose thermometers with a strengthened bulb. A plastic clip prevents rolling.



| Product No                 | Range        | Divisions | Back colour |
|----------------------------|--------------|-----------|-------------|
| 305mm long, 76mm immersion |              |           |             |
| 10452/01                   | -10 to 50°C  | 0.5°C     | White       |
| 10452/02                   | -10 to 110°C | 0.5°C     | White       |
| 10452/03                   | -10 to 250°C | 1.0°C     | White       |
| 10452/04                   | -10 to 360°C | 2.0°C     | White       |
| 10452/05                   | 0 to 240°F   | 1.0°F     | White       |
| 10452/21                   | -10 to 50°C  | 0.5°C     | Yellow      |
| 10452/22                   | -10 to 110°C | 0.5°C     | Yellow      |
| 10452/23                   | -10 to 250°C | 1.0°C     | Yellow      |
| 10452/24                   | -10 to 360°C | 2.0°C     | Yellow      |
| 10452/25                   | 0 to 240°F   | 2.0°F     | Yellow      |

### Spirit Filled

Mercury-free thermometers with a strengthened bulb. A plastic clip prevents rolling

|                             |              |       |        |
|-----------------------------|--------------|-------|--------|
| 152mm long, total immersion |              |       |        |
| 10455/01                    | -10 to 50°C  | 1.0°C | Yellow |
| 10455/02                    | -10 to 110°C | 1.0°C | Yellow |
| 10455/03                    | -10 to 150°C | 1.0°C | Yellow |
| 10455/04                    | -10 to 250°C | 5.0°C | Yellow |
| 10455/05                    | -10 to 360°C | 5.0°C | Yellow |

| Product No | Range          | Divisions | Spirit colour | Immersion length | Total length |
|------------|----------------|-----------|---------------|------------------|--------------|
| 10453/02   | -10 to 110°C   | 1.0°C     | Red           | Total            | 204mm        |
| 10453/03   | -10 to 50°C    | 0.5°C     | Red           | 76mm             | 305mm        |
| 10453/04   | -10 to 110°C   | 1.0°C     | Red           | 76mm             | 305mm        |
| 10453/05   | -10 to 150°C   | 1.0°C     | Red           | 76mm             | 305mm        |
| 10453/06   | -10 to 110°C/F | 1.0°C/F   | Red           | 76mm             | 305mm        |
| 10453/07   | 30 to 220°F    | 1.0°F     | Red           | 76mm             | 305mm        |
| 10453/21   | -10 to 110°C   | 1.0°C     | Green         | Total            | 152mm        |
| 10453/22   | -10 to 110°C   | 1.0°C     | Green         | Total            | 203mm        |
| 10453/23   | -10 to 110°C   | 1.0°C     | Green         | 76mm             | 305mm        |
| 10453/24   | -10 to 110°C/F | 1.0°C/F   | Green         | 76mm             | 305mm        |
| 10453/25   | 0 to 240°F     | 1.0°F     | Green         | 76mm             | 305mm        |
| 10453/41   | -10 to 110°C   | 1.0°C     | Blue          | Total            | 152mm        |
| 10453/42   | -10 to 110°C   | 1.0°C     | Blue          | 76mm             | 305mm        |

BS593, N.A.M.A.S Certified

Supplied complete with U.K.A.S. certificate

| Product No | BS schedule mark | Range         | Divisions | Immersion length | Total length | Ice point scale |
|------------|------------------|---------------|-----------|------------------|--------------|-----------------|
| 10878/01   | A40C/100         | 0 to +40°C    | 0.1°C     | 100mm            | 405mm        | No              |
| 10878/02   | B60C/100         | -20 to +60°C  | 0.2°C     | 100mm            | 405mm        | No              |
| 10878/03   | C105C/50         | -5 to +105°C  | 1°C       | 50mm             | 205mm        | No              |
| 10878/04   | C105C/100        | -5 to +105°C  | 1°C       | 100mm            | 255mm        | No              |
| 10878/05   | F50C/100         | -2 to +52°C   | 0.2°C     | 100mm            | 385mm        | No              |
| 10878/06   | F150C/100        | +98 to +152°C | 0.2°C     | 100mm            | 385mm        | No              |
| 10878/07   | F300C/76         | -5 to +300°C  | 1°C       | 76mm             | 383mm        | No              |
| 10878/08   | A40C/Total       | 0 to +40°C    | 0.1°C     | Total            | 405mm        | No              |
| 10878/09   | B60C/Total       | -20 to +60°C  | 0.2°C     | Total            | 405mm        | No              |