


# Thermometers and Timers

## Temperature Indicating Strips

### Irreversible

A low-cost and simple means of checking the maximum temperature reached, such as in sterilisation procedures. Each self-adhesive and waterproof strip is divided into five or eight temperature steps. As the temperature rises each step (or segment) turns irreversibly black to provide a permanent record. Both Centigrade and Fahrenheit scales are shown and temperatures are accurate  $\pm 2\%$ . Supplied in packs of 10



### Product No Temperature Steps (°C)

11443/01 40, 43, 46, 49, 54

11443/02 60, 65, 71, 77, 82

11443/03 88, 93, 99, 104, 110

11443/04 116, 121, 127, 132, 138

11443/05 143, 149, 154, 160, 166

11443/06 171, 177, 182, 188, 193

11443/07 199, 204, 210, 216, 224

11443/08 232, 241, 249, 254, 260

11443/09 40, 43, 46, 49, 54, 60, 65, 71

11443/10 77, 82, 88, 93, 99, 104, 110, 116

11443/11 121, 127, 132, 138, 143, 149, 154, 160

11443/12 166, 171, 177, 182, 188, 193, 199, 204

11443/13 210, 216, 224, 232, 241, 249, 254, 260

We can also supply single step strips at any one of the above temperatures

## Timers



### Model no TM-22

- 4 (T1/T2/T3/T4)timers in one unit
- 99:59:59
- Count up and count down function

Product no 14010/01



### Model no TM-44

- 4 (T1/T2/T3/T4)timers in one unit
- 99:59:59
- Count up and count down function
- Alarm clock with snooze setting
- 12/24 selectable

Product no 14010/02



### Model no TM-100

- 3 (T1/T2/T3/T4)timers in one unit
- 23:59:59
- Triple row display for easy reading
- Count up and count down function
- Alarm clock with snooze setting
- 12/24 selectable

Product no 14010/03



### Model no TM-110

- 3 (T1/T2/T3/T4)timers in one unit
- 23:59:59
- Triple row display for easy reading
- Count up and count down function
- Alarm clock with snooze setting
- 12/24 selectable

Product no 14010/04



### Model no TM-120

- 3 (T1/T2/T3/T4)timers in one unit
- 23:59:59
- Triple row display for easy reading
- Count up and count down function
- Alarm clock with snooze setting
- 12/24 selectable

Product no 14010/05