

Bunsen Tubing

Black neoprene tubes with reinforced ends for repeated connection/disconnection of methane burners. All lengths have a 7-8mm bore

Product No	Length
10340/01	450mm
10340/02	600mm
10340/03	750mm
10340/04	900mm

Peristaltic Tubing

PVC

- Withstands temperatures from -60 to +70°C
- Food grade (FDA approved)
- Offers good clarity and smooth bores
- Wall thickness is 1.6mm
- Supplied in coils of 15 metres

Product No	Bore	Wall	O. Diameter
11413/01	1.60mm	1.6mm	4.80mm
11413/02	2.40mm	1.6mm	5.60mm
11413/03	3.20mm	1.6mm	6.40mm
11413/04	4.80mm	1.6mm	8.00mm
11413/05	6.40mm	1.6mm	9.60mm
11413/06	8.00mm	1.6mm	11.20mm
11413/07	9.60mm	1.6mm	12.80mm

Silicone

- Withstands temperatures from -60 to +200°C (up to 220°C for short periods)
- Food grade
- Ozone resistant
- High electrical resistance
- Suitable for steam or dry sterilisation
- Contains no toxic or leachable plasticisers
- High strength for long life
- Supplied in coils of 15 metres

Product No	Bore	Wall	O. Diameter
11414/01	0.8mm	1.6mm	4.0mm
11414/02	1.2mm	1.5mm	4.2mm
11414/03	1.6mm	1.6mm	4.8mm
11414/04	3.2mm	1.6mm	6.4mm
11414/05	4.8mm	1.6mm	8.0mm
11414/06	6.4mm	1.6mm	9.6mm
11414/07	8.0mm	1.6mm	11.2mm
11414/08	9.6mm	1.6mm	12.8mm

Neoprene Tubing

Black

Superior to the Dark Grey grade with improved resistance to kinking and bursting. Ideal for use with oils, greases and solvents, and will withstand temperatures from -30°C to 100°C.

Normal Wall

Product No	Bore	Wall	O.Dia	Coil length
11368/01	3.0mm	0.75mm	4.5mm	10m
11368/02	5.0mm	1.5mm	8.0mm	10m
11368/04	6.5mm	1.5mm	9.5mm	10m
11368/05	8.0mm	2.0mm	12.0mm	5m
11368/06	10.0mm	2.0mm	14.0mm	5m
11368/07	12.5mm	2.25mm	17.0mm	5m
11368/08	16.0mm	3.25mm	22.5mm	5m
11368/09	20.0mm	3.25mm	26.5mm	5m
11368/10	25.0mm	3.25mm	31.5mm	5m

Heavy Wall

Suitable to above, but more suited to vacuum or pressurised applications up to one bar

Product No	Bore	Wall	O.Dia	Coil length
11369/01	3.0mm	3.0mm	9mm	5 m
11369/02	5.0mm	4.0mm	13mm	5 m
11369/03	6.5mm	5.0mm	16.5mm	5 m
11369/04	8.0mm	6.0mm	20.0mm	5 m
11369/05	10mm	7.0mm	24.0mm	5 m
11369/06	12.5mm	8.5mm	29.5mm	5 m

PTFE Tubing

- Withstands temperatures from -200 to +200°C
- Food grade
- Good radiation resistance
- Good ozone resistance
- Good UV resistance
- High electrical resistance
- Excellent chemical resistance
- Non-flammable
- Suitable for steam or dry sterilisation
- Supplied in 15 metre coils

Product No	Bore	Wall	O.Dia
11416/01	0.97mm	0.3mm	1.57mm
11416/02	1.01mm	0.3mm	1.61mm
11416/03	1.19mm	0.3mm	1.79mm
11416/04	1.35mm	0.3mm	1.95mm
11416/05	1.50mm	0.3mm	2.10mm
11416/06	1.68mm	0.3mm	2.28mm
11416/07	1.93mm	0.3mm	2.53mm
11416/08	2.16mm	0.3mm	2.76mm
11416/09	2.41mm	0.3mm	3.01mm
11416/10	2.69mm	0.3mm	3.29mm