

# Tubing

## PVC Tubing

Very flexible, non-toxic and resistant to temperatures from -60° to 80°C



Product No	Bore	Wall	Coil length
10283/01	3.0mm	0.75mm	30m
10283/02	5.0mm	1.50mm	30m
10283/03	6.5mm	1.50mm	30m
10283/04	8.0mm	2.00mm	30m
10283/05	10.0mm	2.00mm	30m
10283/06	12.5mm	2.25mm	10m
10283/07	16.0mm	3.25mm	10m
10283/08	20.0mm	3.25mm	10m
10283/09	25.0mm	3.25mm	10m

### Braided

PVC tubing with terylene braiding giving superior strength, particularly under pressure



Product No	Bore	Wall	Coil length
10896/01	5.0mm	3mm	30m
10896/02	6.5mm	3mm	30m
10896/03	8.0mm	3mm	30m
10896/04	10.0mm	3mm	30m
10896/05	12.5mm	3mm	30m
10896/06	16.0mm	3mm	30m
10896/07	19.0mm	4mm	30m

### Medical Grade

High quality, non-toxic, flexible to BS Softness 52 and kink-resistant. Tubing may be autoclaved providing standard cycles are used and care is taken to protect the tubing from stress while hot. Working pressure ranges from 14 to 39psi, dependent upon tubing size. Supplied in coils of 30 metres

### Translucent

Product No	Bore	Wall	O. Diameter
11419/01	0.63mm	0.385mm	1.40mm
11419/02	1.00mm	0.50mm	2.00mm
11419/03	2.00mm	0.50mm	3.00mm
11419/04	3.00mm	0.625mm	4.25mm
11419/05	3.50mm	0.75mm	5.00mm
11419/06	4.00mm	1.00mm	6.00mm
11419/07	5.00mm	1.00mm	7.00mm
11419/08	6.00mm	0.80mm	7.60mm
11419/09	6.50mm	3.10mm	12.70mm
11419/10	7.00mm	1.25mm	9.50mm
11419/11	8.00mm	1.50mm	11.00mm
11419/12	9.00mm	1.60mm	12.20mm

## PVC Tubing

### Clear

Product No	Bore	Wall	O. Diameter
11418/01	1.00mm	0.30mm	1.60mm
11418/02	1.50mm	0.30mm	2.10mm
11418/03	2.00mm	0.50mm	3.00mm
11418/04	2.50mm	1.90mm	6.30mm
11418/05	3.00mm	0.60mm	4.20mm
11418/06	3.50mm	1.50mm	6.50mm
11418/07	4.00mm	1.40mm	6.80mm
11418/08	4.00mm	0.75mm	5.50mm
11418/09	4.50mm	0.90mm	6.30mm
11418/10	5.00mm	1.00mm	7.00mm
11418/11	6.00mm	0.75mm	7.50mm
11418/12	6.35mm	1.575mm	9.50mm
11418/13	6.60mm	1.45mm	9.50mm
11418/14	8.00mm	1.75mm	11.50mm
11418/15	9.50mm	2.00mm	13.50mm

## Polythene Tubing

### Fine Bore, Medical Grade

Supplied in 30 metre coils

Product No	Bore	Wall	O. Diameter
11415/01	0.28mm	0.61mm	1.50mm
11415/02	0.38mm	1.09mm	2.56mm
11415/03	0.40mm	0.80mm	2.00mm
11415/04	0.50mm	1.00mm	2.50mm
11415/05	0.58mm	0.96mm	2.50mm
11415/06	0.76mm	1.22mm	3.20mm
11415/07	0.86mm	1.27mm	3.40mm
11415/08	0.86mm	1.52mm	3.90mm
11415/09	1.02mm	1.98mm	4.98mm
11415/10	1.14mm	1.57mm	4.28mm
11415/11	1.19mm	1.70mm	4.59mm
11415/12	1.40mm	1.90mm	5.20mm
11415/13	1.50mm	2.70mm	6.90mm
11415/14	1.57mm	2.08mm	6.73mm
11415/15	1.67mm	2.42mm	6.51mm
11415/16	2.00mm	3.00mm	8.00mm
11415/17	2.16mm	3.25mm	8.66mm
11415/18	2.45mm	3.70mm	9.85mm
11415/19	2.80mm	4.10mm	11.00mm
11415/20	3.00mm	3.88mm	10.76mm
11415/21	4.00mm	6.00mm	16.00mm