

## Rubber Tubing

Natural red rubber tubing resistant to temperatures from -40°C to 150°C. Supplied in coils of 10 metres



### Normal Wall

Product No	Bore	Wall	O. Diameter
10280/01	3.0mm	0.75mm	4.5mm
10280/02	5.0mm	1.5mm	8.0mm
10280/03	6.5mm	1.5mm	9.5mm
10280/04	8.0mm	2.0mm	12.0mm
10280/05	10mm	2.0mm	14.0mm
10280/06	12.5mm	2.25mm	17.0mm
10280/07	16mm	3.25mm	22.5mm
10280/08	20mm	3.25mm	26.5mm
10280/09	25mm	3.25mm	31.5mm

### Heavy Wall

Similar to above but more suited to vacuum or pressurised applications. Supplied in coils of 10 metres

Product No	Bore	Wall	O. Diameter
10281/01	3.0mm	3.0mm	9.0mm
10281/02	5.0mm	4.0mm	13.0mm
10281/03	6.5mm	5.0mm	16.5mm
10281/04	8.0mm	6.0mm	20.0mm
10281/05	10mm	7.0mm	24.0mm
10281/06	12.5mm	8.5mm	29.5mm

## Viton Tubing

Opaque, black tubing, offering excellent resistance to corrosives, solvents and oils at elevated temperature and with low gas permeability. Suitable for use at temperatures from -31 to 205°C. May be sterilised at 249°C. Supplied in lengths of 7.6 metres.

Product No	Nominal bore	O. Dia	Wall
11420/01	0.8mm	4.1mm	1.6
11420/02	1.6mm	4.9mm	1.6
11420/03	3.1mm	6.4mm	1.6
11420/04	4.8mm	8.3mm	1.7
11420/05	4.8mm	9.9mm	2.5
11420/06	6.4mm	9.9mm	1.7
11420/07	6.4mm	11.2mm	2.4
11420/08	7.9mm	11.2mm	1.7
11420/09	7.9mm	-	-
11420/10	9.7mm	-	-

## Silicone Tubing

### Medical Grade

Non-toxic silicone tubing, approved by the National Water Council. Supplied by the metre with uncut lengths available up to 10 metres if requested. Silicone will withstand repeated sterilisation.

Product No	Bore	Wall
10328/01	1.0mm	1.0mm
10328/02	1.5mm	0.5mm
10328/03	2.0mm	1.0mm
10328/04	3.0mm	0.75mm
10328/05	3.0mm	1.0mm
10328/06	3.2mm	0.8mm
10328/07	3.2mm	1.6mm
10328/08	4.0mm	1.0mm
10328/09	4.5mm	1.6mm
10328/10	5.0mm	1.5mm
10328/11	6.0mm	1.5mm
10328/12	6.3mm	1.6mm
10328/13	6.3mm	3.2mm
10328/14	6.5mm	1.5mm
10328/15	7.9mm	3.2mm
10328/16	8.0mm	1.6mm
10328/17	8.0mm	2.0mm
10328/18	9.0mm	2.0mm
10328/19	9.5mm	3.2mm
10328/20	10.0mm	2.0mm
10328/21	12.5mm	2.25mm
10328/22	12.5mm	3.25mm
10328/23	16.0mm	3.25mm
10328/24	20.0mm	3.25mm
10328/25	25.0mm	3.25mm

