


Tubing Accessories

Tubing Connectors

All (except the Disconnectors) are made of polypropylene, having excellent chemical and temperature resistance.



Straight, Tapered

Tapered to fit a range of tubing bores. Supplied in packs of 10.

Product No	Nominal O. Dia	Internal Dia
10039/01	4/5/6mm	2mm
10039/02	6/7/8mm	3.5mm
10039/03	8/9/10mm	4.5mm
10039/04	10/11/12mm	6.5mm
10039/05	12/13/14mm	8.5mm
10039/06	14/15/16mm	10.5mm

'T' Pattern

Ridges help prevent slippage and aid sealing.



Product No	Nominal O. Dia	Internal Dia
10417/01	4mm	1.6mm
10417/02	6mm	2.8mm
10417/03	8mm	4.4mm
10417/04	10mm	6.3mm
10417/05	12mm	8.2mm
10417/06	14mm	10.0mm
10417/07	16mm	12.1mm

'Y' Pattern



Product No	Nominal O. Dia	Internal Dia
10327/01	4mm	2.0mm
10327/02	6mm	2.7mm
10327/03	8mm	4.6mm
10327/04	10mm	5.5mm
10327/05	12mm	7.3mm
10327/06	14mm	9.7mm
10327/07	16mm	11.9mm

Tubing Connectors

Unequal, Tapered

For connecting tubing of differing Internal Diameters. Outer and Internal Diameters are given for the small & large connections, respectively.



Product No	Nominal O. Dia	Internal Dia
10420/01	4/6/8 & 10/12mm	1.6 & 4.6mm
10420/02	4/6/8 & 12/14/16mm	1.6 & 7.5mm
10420/03	8/10/12 & 14/16mm	7.1/12.0 & 11.3/15.9mm

Disconnectors, Tapered

Made of polyethylene, these allow quick and simple connection and disconnection without causing undue wear to tubing ends. Tapered to fit a range of tubing bores.



Product No	Nominal O. Dia	Internal Dia
10419/01	4/5/6mm	1.6mm
10419/02	6/7/8mm	2.7mm
10419/03	8/9/10mm	3.6mm
10419/04	10/11/12mm	5.5mm
10419/05	12/13/14mm	7.3mm
10419/06	14/15/16mm	8.8mm

Connectors

CAP / BARB

Straight or angled polypropylene connectors with silicone rubber seals. GL14 aperture caps will be required and should be ordered separately. Supplied in packs of 10.

Product No	Description	I. Dia	O. Dia
10632/01	Straight connectors	5mm	9mm
10632/02	Angled connectors	4mm	8mm
10090/01	Aperture caps GL14	-	-