

Tubing Accessories

Connectors

BARB / BARB

Double ended connectors
115mm long of borosilicate
glass, ideal for viewing
contents of opaque tubing



Product No	Outer Dia	Internal Dia
10630/01	6mm	3mm
10630/02	8mm	5mm
10630/03	9mm	6mm
10630/04	10mm	7mm
10630/05	11mm	8mm
10630/06	12mm	9mm

Non-Return Valves

Will prevent back flow of liquid
in tubing with low vacuum or
up to 2 bar pressure (at 20°C).
Tapered connectors to fit a
range of tubing bores.




Product No	Nominal O. Dia	Internal Dia
10418/01	8/9/10mm	3mm
10418/02	12/13/15mm	5.5mm

Flow Indicators

Suitable for use with non-corrosive liquids and gases up to
85°C. Tapered connectors are ridged to aid tube grip and
sealing

Ball Type


Product No	11421/01
Min. liquid flow	50ml/minute (horizontal) 140ml/minute (vertical)
Min. gas flow	3 litres/minute (horizontal) 8 litres/minute (vertical)
Tube connectors	6.5-10mm Outer Diameter
Overall length	102mm



Flow Indicators

Paddle Wheel Type


Product No	11421/02
Flow range	5 to 95ml/second (0.3 to 5.7 litres/minute)
Max. pressure	2 bar
Tube connectors	6.0-11mm Outer Diameter
Overall length	90mm



Stopcocks

Polypropylene & HDPE Key


Operate smoothly requiring no
lubrication. Suitable for liquids and
gases at up to 0.5 bar and low vacuum



Product No	Nominal Outer Dia	Internal Dia
10421/01	8mm	4.9mm
10421/02	10mm	7.0mm
10421/03	13mm	10.9mm
10421/04	16mm	12.7mm

Polypropylene & PTFE Key

Tapered connectors are ridged for
improved grip and to aid sealing. Will
accept 1/4" to 5/16" tubing.
Autoclavable



Product No	Style	Bore
11422/01	2-Way	2mm
11422/02	2-Way	4mm
11422/03	3-Way	2mm
11422/04	3-Way	4mm

