

Basins, Crucibles and Dishes

Crucibles and Lids

Product no	Dia.	Material	Description	Capacity	Thickness
12301/01	20mm	Nickel	Crucible	5ml	0.5mm
12301/02	25mm	Nickel	Crucible	10ml	0.5mm
12301/03	30mm	Nickel	Crucible	15ml	0.5mm
12301/04	35mm	Nickel	Crucible	25ml	0.5mm
12301/05	40mm	Nickel	Crucible	30ml	0.5mm
12301/06	45mm	Nickel	Crucible	50ml	0.5mm
12301/07	50mm	Nickel	Crucible	70ml	0.5mm
12301/08	60mm	Nickel	Crucible	130ml	0.5mm
12301/09	80mm	Nickel	Crucible	270ml	0.5mm
12301/10	25mm	Nickel	Crucible	10ml	1.0mm
12301/11	30mm	Nickel	Crucible	15ml	1.0mm
12301/12	35mm	Nickel	Crucible	25ml	1.0mm
12301/13	40mm	Nickel	Crucible	30ml	1.0mm
12301/14	45mm	Nickel	Crucible	50ml	1.0mm
12301/15	50mm	Nickel	Crucible	70ml	1.0mm
12301/16	60mm	Nickel	Crucible	130ml	1.0mm
12301/17	80mm	Nickel	Crucible	270ml	1.0mm
12301/18	100mm	Nickel	Crucible	600ml	1.0mm
12301/19	30mm	Nickel	Crucible	15ml	1.0mm
12301/20	40mm	Nickel	Crucible	30ml	1.5mm
12301/21	45mm	Nickel	Crucible	50ml	1.5mm
12301/22	40mm	Nickel	Crucible	30ml	2.0mm
12301/23	45mm	Nickel	Crucible	50ml	2.0mm
12301/24	50mm	Nickel	Crucible	70ml	2.0mm
12301/25	60mm	Nickel	Crucible	130ml	2.0mm
12301/30	20mm	Nickel	Crucible Lid		
12301/31	25mm	Nickel	Crucible Lid		
12301/32	30mm	Nickel	Crucible Lid		
12301/33	35mm	Nickel	Crucible Lid		
12301/34	40mm	Nickel	Crucible Lid		
12301/35	45mm	Nickel	Crucible Lid		
12301/36	50mm	Nickel	Crucible Lid		
12301/37	60mm	Nickel	Crucible Lid		
12301/38	80mm	Nickel	Crucible Lid		
12301/39	100mm	Nickel	Crucible Lid		
12301/50	25mm	Zirconium	Crucible	10ml	
12301/51	30mm	Zirconium	Crucible	15ml	
12301/52	35mm	Zirconium	Crucible	25ml	
12301/53	38mm	Zirconium	Crucible	30ml	
12301/54	40mm	Zirconium	Crucible	35ml	
12301/55	45mm	Zirconium	Crucible	40ml	
12301/56	45mm	Zirconium	Crucible	45ml	
12301/57	45mm	Zirconium	Crucible	50ml	
12301/58	45mm	Zirconium	Crucible	60ml	
12301/59	48mm	Zirconium	Crucible	70ml	
12301/70	25mm	Zirconium	Crucible Lid		
12301/71	30mm	Zirconium	Crucible Lid		
12301/72	35mm	Zirconium	Crucible Lid		
12301/73	38mm	Zirconium	Crucible Lid		
12301/74	40mm	Zirconium	Crucible Lid		
12301/75	45mm	Zirconium	Crucible Lid		
12301/76	45mm	Zirconium	Crucible Lid		
12301/77	45mm	Zirconium	Crucible Lid		
12301/78	45mm	Zirconium	Crucible Lid		
12301/79	48mm	Zirconium	Crucible Lid		



Evaporating Basins

Made from porcelain, with round bottom and spout.



Product No	Capacity	Dia x Height-mm	Form
10053/01	35	63 x 25	Shallow
10053/02	50	70 x 28	Shallow
10053/03	75	80 x 32	Shallow
10053/04	95	85 x 35	Shallow
10053/05	170	105 x 42	Shallow
10053/06	285	125 x 50	Shallow
10053/07	420	150 x 55	Shallow
10053/08	1000	200 x 80	Shallow
10053/09	1225	210 x 85	Shallow
10054/01	50	70 x 34	Deep
10054/02	150	97 x 50	Deep
10054/03	250	110 x 52	Deep
10054/04	700	150 x 68	Deep
10054/05	900	167 x 75	Deep

Silica Basins

Translucent glazed quartz basins with spout will withstand continuous use at 1050°C (or 1300°C for short periods).



Product No	Capacity-ml	Dia x Height-mm	Form
10439/01	20	55 x 23	Low
10439/02	45	75 x 27	Low
10439/03	70	87 x 33	Low
10439/04	100	100 x 38	Low
10439/05	200	113 x 47	Low
10440/01	90	80 x 45	Deep
10441/01	25	79 x 15	Flat
10441/02	75	101 x 20	Flat
10441/03	125	131 x 23	Flat