

Reagents

For use with Comparators or PCcheckit Photometers Lovibond offer a large range of reagents for measuring numerous parameters and ranges. We list here just the most popular, so please ask if you would like information on reagents for other tests. NB If two or three reagents are listed for any test, ALL will be required.

Product No	Test	For discs/PCcheckit	Reagent	Pack Qty
10928/24	Aluminium	3/127, PCcheckit	Aluminium No 1 tablets	250
10928/25			Aluminium No 2 tablets	250
10928/01	Ammonia	3/112, 3/113, 3/125	Ammonia No 1 tablets	250
10928/02		PCcheckit	Ammonia No 2 tablets	250
10928/03	Bromine*	3/53A, 3/53B	DPD No 1 tablets	250
10928/26		PCcheckit	DPD Nitrite tablets	250
10928/27	Chloride	PCcheckit	Chloride T1 tablets	250
10928/28			Chloride T2 tablets	250
10928/03	Chlorine Dioxide, Free	3/40ED	DPD No 1 tablets	250
10928/04	Chlorine, Free*	3/40A, E, K, PCcheckit	DPD No 1 tablets	250
10928/05	Chlorine Dioxide, Total	3/40ED	DPD No 3 tablets	250
10928/05	Chlorine, Total*	3/40A, E, K, PCcheckit	DPD No 3 tablets	250
10928/29	Chlorine, Free*	PCcheckit	Liquid DPD 3 x 15ml each solutions A & B	300 tests
10928/30	Chlorine, Free & Total*	PCcheckit	Liquid DPD 2 x 15ml each solutions A, B & C	200 tests
10928/31	Chlorine Dioxide, Total	3/157	Potassium Iodide tablets	250
10928/31	Chlorine, Total*	3/2ARP, 3/2IOD	Potassium Iodide tablets	250
10928/32		PCcheckit	Acidifying GP tablets	250
10928/34	Copper, Low Range	3/106	Copper/zinc LR tablets	250
10928/35	Copper, High Range	3/110	Copper/zinc HR tablets	250
10928/36	Copper	PCcheckit	Copper No 1 tablets	250
10928/37			Copper No 2 tablets	250
10928/38	Hydrazine	3/126, 3/135, PCcheckit	Hydrazine test powder	30 tests
10928/39	Iron, Low Range	3/116, PCcheckit	Iron LR tablets	250
10928/40	Iron, High Range	3/117	Iron HR tablets	250
10928/41	Manganese	3/134	Manganese No 1 tablets	250
10928/42			Manganese No 2 tablets	250
10928/43	Molybdate	3/137	Molybdate No 1 tablets	250
10928/44			Molybdate No 2 tablets	250
10928/45	Nitrate, Low Range	3/124	Nitrate test powder	100 tests
10928/18			Nitrate test tablets	100
10928/46			Nitrite LR tablets	250
10928/20	Nitrate, High Range	3/142	Nitrate No 1 tablets	250
10928/21			Nitrate No.2 tablets	250
10928/46	Nitrite, Low Range	3/103, PCcheckit	Nitrite LR tablets	250
10928/48	Ozone	3/148, PCcheckit	Ozone tablets	250
10928/49	pH	2/1G	Bromocresol purple tablets	250
10928/50	pH	2/1H	Bromothymol blue tablets	250
10928/51	pH	2/1J	Phenol red tablets	250
10928/52	pH	PCcheckit	Phenol red photometer tablets	250
10928/53	pH	2/1L	Thymol blue tablets	250
10928/54	Phosphate, Low Range	3/133, PCcheckit	Phosphate LR No 1 tablets	250
10928/55			Phosphate LR No 2 tablets	250
10928/56	Phosphate, High Range	3/136, 3/70	Phosphate HR tablets	250
10928/57	Phosphate, High Range	PCcheckit	Phosphate P1 tablets	250
10928/58			Phosphate P2 tablets	250
10928/59	Sulphide	3/128	Sulphide No 1 tablets	100
10928/60			Sulphide No 2 tablets	100
10928/61	Silica	3/139	Silica No 1 tablets	250
10928/62			Silica No 2 tablets	250
10928/34	Zinc, Low Range	PCcheckit	Copper/zinc LR tablets	250
10928/64			EDTA tablets	250

* DPD No 1 tablets are used for determining Active Bromine. Combined and Free Bromine can be measured with the addition of DPD Nitrite tablets.

DPD No 1 tablets alone are used for determining Free Chlorine. Both DPD No 1 and 3 tablets are used for testing Free, Combined and Total Chlorine.