


Water and Wire Holders



Other Reagents & Pack Sizes

We can supply all Lovibond reagents. Here are just some of the alternatives available


Product No	Test	Instrument	Description	Pack Qty
10928/07	Chlorine	Nessleriser	DPD No 1 tablets	250
10928/04	Chlorine	Comparator	DPD No 2 tablets	250
10928/22	Chlorine	Nessleriser	DPD No 2 tablets	100
10928/09	Chlorine	Nessleriser	DPD No 3 tablets	250
10928/06	Chlorine	Comparator	DPD No 4 tablets	250
10928/11	Chlorine	Comparator	Potassium Iodide 0.3g tablets	100 (Foil)
10928/12	Chlorine	Comparator	Potassium Iodide 0.3g tablets	100 (Bottle)
10928/13	Chlorine	Comparator	Acidifying GP tablets	100
10928/17	Resazurin in milk	Comparator	Resazurin tablets	100
10928/14	pH	Comparator	Phenol Red tablets	100
10928/15	pH	Comparator	Bromothymol Blue tablets	100
10928/23	pH	Comparator	Cresol Red tablets	250

Optional Lighting Units

Essential for use in areas without a constant, stable light source and recommended to provide a guarantee of uniform lighting conditions. (Good Northern daylight is the optimum natural lighting for colour matching)

Product No	Description
10921/01	Daylight 2000 Lighting Unit (benchtop)
10922/03	Portable Lighting Unit for Comparator 2000+
10922/02	Charger for Portable Lighting Unit

NB The above Portable Lighting Unit is suitable only for the Comparator 2000+. Those who require a Lighting Unit for the original Comparator 2000 should order Cat. No 10922/01 (17 03 50) plus the above charger. The Daylight 2000 Lighting Unit is suitable for use with either model.



Wire Holder

Product No **11439/01**

An aluminium handle with screw-fastening chuck, for platinum wire, disposable loops and needles up to 0.56mm (24swg).
Length 175mm