

# Baisins, Crucibles and Dishes

## Capsules, Silica

Translucent glazed silica circular or rectangular capsules. Will withstand up to 1050° continuously or 1300°C for short periods.



Product No	Capacity-ml	Dia x Height-mm	L x W x D-mm
10442/01	10	40 x 15	
10442/02	20	48 x 15	
10442/03	25	54 x 15	
10442/04	35	69 x 15	
10442/05	45	56 x 28	
10442/06	50	78 x 18	
10442/07	80	75 x 28	
10442/08	100	87 x 28	
10442/09	23	53 x 13	
10443/01	10		60 x 28 x 12
10443/02	20		68 x 39 x 12
10443/03	25		52 x 42 x 18

## Crucibles, Porcelain

Medium and low form crucibles supplied WITHOUT lids which must be ordered separately (see Product. No.10450/xx below).



### Bases

Product No	Capacity-ml	Dia x Height-mm	Form
10449/01	6	30 x 19	Low
10449/02	11	35 x 22	Low
10449/03	15	40 x 25	Low
10449/04	25	45 x 28	Low
10449/05	50	60 x 38	Low
10449/06	85	70 x 44	Low
10449/21	12	35 x 28	Tall
10449/22	20	40 x 32	Tall
10449/23	30	45 x 36	Tall
10449/24	45	50 x 40	Tall
10449/25	80	60 x 48	Tall
10449/26	120	70 x 56	Tall

### Lids



Product No	To fit Crucible	Lid Dia-mm
10450/01	10449/01	30
10450/02	10449/02, 10449/21	35
10450/03	10449/03, 10449/22	40
10450/04	10449/04, 10449/23	45
10450/05	10449/24	50
10450/06	10449/05, 10449/25	60
10450/07	10449/06, 10449/26	70

## Crucibles, Silica

### Bases

Low, tall & broad form translucent glazed crucibles, supplied WITHOUT lids which must be ordered separately (see Cat. No's 10445/xx and 10447/xx, below). Will withstand up to 1050°C continuously or 1300°C for short periods.



Product No	Capacity -ml	Dia x Height-mm	Form
10444/01	5	29 x 19	Low
10444/02	10	38 x 19	Low
10444/03	15	41 x 25	Low
10444/04	25	47 x 28	Low
10444/05	50	57 x 37	Low
10444/06	80	67 x 45	Low
10444/07	150	82 x 51	Low
10446/01	10	27 x 32	Tall
10446/02	20	35 x 38	Tall
10446/03	30	43 x 39	Tall
10446/04	50	51 x 51	Tall
10446/05	200	76 x 81	Tall
10448/01	15	41 x 24	Broad
10448/02	25	48 x 27	Broad
10448/03	40	57 x 37	Broad

### Lids

Choose according to crucible diameter and form. Low form lids include a ring handle.

Product No	Crucible Dia-mm	Crucible Form
10445/01	29	Low
10445/02	41	Low
10445/03	47	Low
10445/04	57	Low
10445/05	67	Low
10445/06	82	Low
10447/01	27	Tall
10447/02	35	Tall
10447/03	43	Tall
10447/04	51	Tall
10447/05	76	Tall