


# Basins, Crucibles and Dishes

## Crucible Tongs

Available in stainless or mild steel, bowed or straight. Larger sizes available to order.



Product No	Length-mm	Material	Style
10323/01	150	S/steel	Bowed
10323/02	200	S/steel	Bowed
10323/03	300	S/steel	Bowed
10324/01	150	Mild steel	Bowed
10324/02	200	Mild steel	Bowed
10324/03	300	Mild steel	Bowed
10321/01	150	S/steel	Straight
10321/02	200	S/steel	Straight
10321/03	300	S/steel	Straight
10322/01	150	Mild steel	Straight
10322/02	200	Mild steel	Straight
10322/03	300	Mild steel	Straight

## Furnace Tongs

Product no	Description	Length
14020/01	Mild steel, straight with flat hinge	380mm
14020/02	Mild steel, with bow and flat hinge	380mm



## Furnace Tongs, ward pattern



Stainless steel, for handling a wide variety of crucibles, etc, especially in muffle furnaces. Operated by thumb action with one hand and fitted with a circular shield to protect the hand against radiation from the furnace. Dimensions given are overall length x length to shield

Product no	Length	Length to shield
14020/03	590mm	410mm
14020/04	740mm	560mm

## Flask & Beaker Tongs

**Beaker** Stainless steel, with 25mm wide jaws lined with replaceable glass fibre which reduces the risk of breakage. Firmly grips beakers from 100 to 1,500mL.



Product no	Length mm
14020/10	325

### Spares

Product no	Description	Pack qty
14020/12	Glass fibre sleeves for 14020/10	2

**Flask** Stainless steel to hold flasks from 100 to 2000mL.

Product no	Length mm
14020/11	230