

# Autoclaves and Accessories

## Autoclave Baskets



### Baskets, Polypropylene

- Autoclavable at 121°C
- Ideal for sterilising test tubes
- Useful for storage in the laboratory
- Stackable, space saving design

| Product no      | l x w x h mm    | Material |
|-----------------|-----------------|----------|
| <b>12008/01</b> | 143 x 120 x 109 | PP       |
| <b>12008/02</b> | 162 x 162 x 155 | PP       |
| <b>12008/03</b> | 225 x 225 x 230 | PP       |

## Autoclave Deodouriser Capsules

Gelatine capsules containing highly concentrated deodourising materials, specially formulated to neutralise the pungent malodours which can occur when autoclaving equipment is soiled with organic matter. Suitable for use in all makes of laboratory, medical and industrial autoclaves.



| Product no      | Pack Qty |
|-----------------|----------|
| <b>12042/20</b> | 100      |
| <b>12042/21</b> | 500      |

## Autoclave Baskets

White nylon-coated mild steel wire.  
(Also available in electro-polished stainless steel).

| Product no      | Cube Size |
|-----------------|-----------|
| <b>10055/01</b> | 100mm     |
| <b>10055/02</b> | 125mm     |
| <b>10055/03</b> | 150mm     |
| <b>10055/04</b> | 200mm     |

## Autoclave Bottle

Wide-mouth, square, polypropylene 4-litre bottle with moulded handgrips, and polypropylene cap. Excellent for storing and sterilising contaminated slides and other small items.



| Product no | 11009/01 |
|------------|----------|
|            |          |

## Bags, Autoclave, 121°C

Specifically designed for contaminated waste disposal in autoclaves or incinerators. Suitable for sterilisation at 121°C and at 135°C for the decontamination and inactivation of particularly resistant biological waste. Supplied in convenient "tissue box" cartons dispense bags individually



| Product No      | Description   | Max Temp | Case      | Qty |
|-----------------|---------------|----------|-----------|-----|
| Width x Length  |               | (°C)     | (mm)      |     |
| <b>12008/20</b> | Autoclave Bag | 121      | 305 x 660 | 200 |
| <b>12008/21</b> | Autoclave Bag | 121      | 305 x 660 | 500 |
| <b>12008/22</b> | Autoclave Bag | 121      | 406 x 610 | 200 |
| <b>12008/23</b> | Autoclave Bag | 121      | 406 x 610 | 500 |
| <b>12008/24</b> | Autoclave Bag | 121      | 610 x 810 | 200 |
| <b>12008/25</b> | Autoclave Bag | 135      | 307 x 660 | 200 |
| <b>12008/26</b> | Autoclave Bag | 135      | 307 x 660 | 500 |
| <b>12008/27</b> | Autoclave Bag | 135      | 406 x 610 | 200 |
| <b>12008/28</b> | Autoclave Bag | 135      | 406 x 610 | 500 |
| <b>12008/29</b> | Autoclave Bag | 135      | 610 x 810 | 200 |

## Bags, Autoclave, Holders

Coated wire or cardboard holders for use with autoclave bags



Cardboard Holder



Coated Wire Holder

| Product no      | Description               | Material    | Case Qty |
|-----------------|---------------------------|-------------|----------|
| <b>12008/40</b> | Holder for 511 bags       | Cardboard   |          |
| <b>12008/41</b> | Holder for 509 & 510 bags | Coated Wire |          |
| <b>12008/42</b> | Holder for 511 bags       | Coated Wire | 1        |

## Autoclave Indicator Tapes

Self-adhesive, with bands that change from white to dark brown when exposed to steam for 3 1/2 minutes at 134°C or 20 minutes at 121°C. Supplied in 55 metre rolls.

| Product no      | Width |
|-----------------|-------|
| <b>10017/01</b> | 19mm  |
| <b>10017/02</b> | 25mm  |