

## Unstirred bath

A range of basic baths with analogue control. A clip is provided to hold a reference thermometer for accurate temperature setting and monitoring. A polypropylene tank, with outer coating of PVC-coated zinc treated steel gives excellent corrosion resistance. Straight forward operating instructions are printed on the front panel. Perforated tray allows mixing of liquid by convection, and for safe operation all except the 10062/01 include an over-temperature cut-out

| Product No | Capacity (Litres) |
|------------|-------------------|
| 10062/01   | 3.5               |
| 10062/02   | 10                |
| 10062/03   | 16                |
| 10062/04   | 24                |
| Product No | Lids to fit       |
| 10063/01   | 10062/01          |
| 10063/02   | 10062/02          |
| 10063/03   | 10062/03          |
| 10063/04   | 10062/04          |



## Transparent Bath

A 3.5 litre bath having a clear polycarbonate tank, with a removable control unit and perforated tray. A user-resettable fixed temperature cut-out is provided for safety. The bath has a temperature range of ambient +5°C to 60°C ±0.3°C.

| Product No | Capacity (Litres) |
|------------|-------------------|
| 10064/01   | 3.5               |
| Product No | Lids to fit       |
| 10065/01   | 10064/01          |



## Shaking Bath

User can select orbital or linear shaking. Operates at ambient to 99°C, or from 0-99°C using optional CS200G Refrigerated Immersion Cooler. Precise control of speed and temperature via digital control and LED display. Racks and trays should be ordered separately.

| Product No          | Description                       |
|---------------------|-----------------------------------|
| 10978/01            | Orbital/Linear Shaking Water Bath |
| Specification       |                                   |
| Tank Size (L, W, D) | 505 x 300 x 200                   |
| Tray Size           | 375 x 235                         |
| Shaking Speed       | 20-200 rpm (orbital)              |
| Shaking Speed       | 40-360 strokes/minute (linear)    |

### Accessories

| Product No | Description                          |
|------------|--------------------------------------|
| 10978/10   | Universal Flask Tray for 10978/01    |
| 10978/11   | Tray to hold five H1 Test Tube Racks |
| 10978/12   | Rack for 48 x 10mm Tubes             |
| 10978/13   | Rack for 44 x 13mm Tubes             |
| 10978/14   | Rack for 24 x 16mm Tubes             |
| 10978/15   | Rack for 21 x 19mm Tubes             |
| 10978/16   | Rack for 12 x 25mm Tubes             |
| 10978/17   | Rack for 10 x 30mm Tubes             |
| 10978/18   | Perforated Tray                      |
| 10978/19   | Plain tray                           |
| 10978/20   | Stainless Steel Gabled Lid           |

