

Ultrasonic Bath

Ideal for cleaning, de-gassing, dissolution and sonochemistry. The tanks are of stainless steel with radiused corners. The XB14 and XB22 models have a variable timer (up to 30 minutes), a separate heater switch to operate typically at 50°C, a resettable safety cut-out operating at 75°C or in case of low liquid level. The XB14 and XB22 have a fitted drainage tap. Baskets and lids must be ordered separately. Flat lids of stainless steel minimise operating noise and evaporation.



Product No	Capacity (Litres)	Working Volume
10977/01	2	149 x 135 x 100
10977/02	3	238 x 135 x 100
10977/03	6	150 x 300 x 140
10977/04	14	325 x 300 x 140
10977/05	22	505 x 300 x 140

Ultrasonic Bath - Accessories

Baskets of nylon coated stainless steel, with lifting handles

Product No	To Fit
10977/10	10977/01
10977/11	10977/02
10977/12	10977/03
10977/13	10977/04
10977/14	10977/05

Flat lids of stainless steel (minimise noise & evaporation)

Product No	To Fit
10977/20	10977/01
10977/21	10977/02
10977/22	10977/03
10977/23	10977/04
10977/24	10977/05

