

Baths, Accessories

Accessories for Water Baths and Shaking Water Baths

- The lift up Makrolon cover is recommended for use at working temperatures above +60 °C for both models
- Large selection of carrier trays for Erlenmeyer flasks and test tubes is available
- For counter-cooling at or near ambient temperature, use the cooling coil (product no **12250/03**)
- To maintain a constant liquid level use the liquid level/cooling set (product no **12250/02**)

Bath Cover

Product code	Description
12250/01	Lift up Makrolon cover

Cooling installations/continuous water supply

Product code	Description
12250/02	Liquid Level/cooling set
12250/03	Cooling coil (figure 1)



Standard carrier trays

Product code	Description
12250/50	For 45 flasks 25 ml
12250/51	For 32 flasks 50ml
12250/52	For 18 flasks 100ml
12250/53	for 15 flasks 200m
12250/54	for 11 flasks 250-300ml
12250/55	for 8 flasks 500ml
12250/56	for 5 flasks 1000ml



All purpose spring tray

Product code	Description
12250/05	pre assembled for 11 Erlenmeyer flasks 250ml including set of springs for Erlenmeyer flasks 25...1000ml



Carrier trays and spring clamps

Product code	Description
12250/06	Basic tray for assembling spring clamps on a mix and match basis

Spring clamps

Product code	Description
12250/10	for 10ml flasks
12250/11	for 25ml flasks
12250/12	for 50ml flasks
12250/13	for 100ml flasks



Product code	Description
12250/60	for 200-250ml flasks
12250/61	for 300ml flasks
12250/62	for 500ml flasks
12250/63	for 1000ml flasks

Picture shows the basic tray assembled with different spring clamps

Carrier tray with test tube racks

12250/20	Basic tray for assembling a maximum of 4 test tubes
-----------------	---

Test tube racks - made of Polypropylene 0 (to +80 °C)

Product code	Description
12250/25	for 60 test tubes, 16/17mm dia.
12250/26	for 90 test tubes, 12/13mm dia.
12250/27	for 90 microliter tubes, 11/12mm dia.
12250/28	for 21 test tubes, 30mm dia.



Test tube racks - made of stainless steel (to +100°C)

Product code	Description
12250/70	for 50 test tubes, 16/17mm dia.
12250/71	for 90 test tubes, 12/13mm dia.
12250/72	for 90 microliter tubes, 11/12mm dia.
12250/73	for 21 test tubes, 30mm dia.

Carrier trays with test tube racks (complete)

(to +80°C)

Product code	Description
12250/30	for 240 test tubes, 16/17mm dia.
12250/31	for 360 test tubes, 12/13mm dia.
12250/32	for 360 microliter tubes, 30 x 11/12mm dia
12250/33	for 84 test tubes, 30mm dia.

Lab automation

Product code	Description
12250/80	'Easy Temp' control software (free download www.julabo.de)
12250/81	RS232 interface cable, 3m for direct PC connection
12250/82	USB interface adapter cable

Option	Description
12250/35	Shaking tray positioner The shaking water bath must be factory-prefitted with the shaking tray positioner. Please specify when placing an order. it cannot be rectified.