

Balances

Trooper™ Basic - Economical compact industrial bench scales

Applications


- Weighing (4 units)
- Mains or battery operation

Features

- RS232 interface
- Backlit LCD display
- Stainless steel platform
- Auto-timer shut-off
- Levelling feet and bubble
- Refer to later pages for the correct calibration weights

Options

- EC Type Approval - 2nd stage



Product no	Capacity	Factory readability	Approved Readability	Platform
12205/60	3kg	0.5g	1g	340 x 230mm
12205/61	6kg	1g	2g	340 x 230mm
12205/62	15kg	2g	5g	340 x 230mm
12205/63	30kg	5g	10g	340 x 230mm

Options - Must be ordered at the same time as the scales

Product no	Description
12205/64	2nd Stage verification* (add 1 to the end of the model no) *postcode of final destination must be supplied at time of order

Accessories

Product no	Description
12205/65	Carrying case
12205/66	PC interface cable - 9 pin
12205/67	Printer SF42 (includes interface cable)
12205/68	Replacement AC mains adaptor

