

Beakers

Beakers, Low Form Manufactured of borosilicate heat resistant glass to ISO3819, DIN12331. With spout and graduations.

Product No	Capacity (ml)	Std pack quantity
10024/01	5	10
10024/02	10	10
10024/03	25	10
10024/04	50	10
10024/05	100	10
10024/15	150	10
10024/06	250	10
10024/07	400	10
10024/08	600	10
10024/09	800	10
10024/10	1000	10
10024/11	2000	10
10024/12	3000	5
10024/13	5000	1
10024/14	10000	1



Beakers, Tall Form Manufactured of borosilicate heat resistant glass to ISO3819, DIN12331. With spout and graduations.

Product No	Capacity (ml)	Std Pack quantity
10070/01	50	10
10070/02	100	10
10070/03	150	10
10070/04	250	10
10070/05	400	10
10070/06	600	10
10070/07	800	10
10070/08	1000	10
10070/09	2000	10
10070/10	3000	2



Beakers, Heavy Duty Heavy duty beakers with spout. Made of borosilicate glass with reinforced rims and thicker walls & bases for increased strength and longer life. With graduations, except the 10, 15 and 20 litre sizes.

Product No	Capacity (ml)	Std Pack quantity
10071/08	100	12
10071/01	150	12
10071/02	250	12
10071/03	400	12
10071/09	500	12
10071/04	600	12
10071/05	1000	12
10071/06	2000	4
10071/10	3000	4
10071/07	4000	4
10071/11	5000	4
10071/12	10000	1
10071/13	15000	1
10071/14	20000	1

