

Beakers, polypropylene

Moulded in one piece from robust polypropylene with spout and graduations in ml and oz. Autoclavable to 121°C for 20 mins.



Product No	Capacity-ml	Pack Qty
10072/01	25	10
10072/02	50	10
10072/03	100	10
10072/04	250	10
10072/05	500	5
10072/06	1000	1
10072/07	2000	1
10072/08	3000	1
10072/09	5000	1

Beakers, Blue Graduations

Manufactured from polypropylene to give excellent chemical resistance and, with permanent blue graduations ($\pm 10\%$) these beakers conform to ISO7056 - 1981 (E) and BS5404, Part 1. They offer very high translucency, and are autoclavable to 121°C for 20 minutes (or continuously at 100°C).



Product No	Capacity x Div-ml	Dia x Height-mm	Pack Qty
11031/01	25 x 1	34 x 49	10
11031/02	50 x 2	41 x 60	10
11031/03	100 x 5	51 x 72	10
11031/04	250 x 10	71 x 95	10
11031/05	500 x 10	87 x 119	5
11031/06	1000 x 20	109 x 147	1
11031/07	2000 x 50	132 x 184	1
11031/08	3000 x 500	158 x 201	1
11031/09	5000 x 500	189 x 229	1

Beakers, Polypropylene, Tripour

A novel design provides three dripless pouring spouts. Supplied in packs of 10 beakers.

- Resistant to most acids, alkalis and solvents
- Autoclavable to 121°C
- Stackable to save storage space
- Graduated ($\pm 5\%$)

Wax-impregnated Tri-pour beaker caps are available on special order in packs of 250 & 500.

Product No	Capacity-ml	Divisions
11157/01	50	5
11157/02	100	10
11157/03	250	10
11157/04	400	20
11157/05	800	50
11157/06	1000	50

Beakers, TPX

Moulded in one piece from heat and chemical resistant TPX with spout and moulded graduations in ml and oz. will withstand 175°C for short periods, 150°C for continuous work.



Product No	Capacity-ml	Pack Qty
10073/01	25	10
10073/02	50	10
10073/03	100	10
10073/04	250	10
10073/05	500	5
10073/06	1000	11
10073/07	2000	1
10073/08	3000	1
10073/09	5000	1

Beakers, TPX, Blue Graduations

As above, but with permanent blue graduations.

Product No	Capacity x Div	OD x Height-mm	Pack Qty
11032/01	25 x 1	34 x 49	10
11032/02	50 x 2	42 x 60	10
11032/03	100 x 5	52 x 72	10
11032/04	250 x 10	71 x 95	10
11032/05	500 x 10	88 x 119	5
11032/06	1000 x 20	110 x 146	1
11032/07	2000 x 50	133 x 184	1
11032/08	3000 x 500	159 x 201	1
11032/09	5000 x 500	190 x 228	1

Beakers, PTFE

Moulded PTFE with a flat base for even heat transfer. Highly chemically resistant, inert, and usable to 280°C. Supplied singly in the sizes below.



Product No	Capacity-ml	Height-mm	Diameter-mm
10733/01	1	18	13
10733/02	5	26	20
10733/03	10	33	24
10733/04	25	47	32
10733/05	50	60	43
10733/06	100	68	54
10733/07	150	69	59
10733/08	250	97	66
10733/09	400	106	80
10733/10	500	125	80
10733/11	600	125	90
10733/12	1000	155	100
10733/13	2000	210	120
10733/14	5000	280	170