

Bottles and Jars

Powder Jars

Clear Uncapped

Clear glass, extra wide mouth

Product No	Capacity-ml	Cap Size
10937/01	15	R3/28
10937/02	30	R3/33
10937/03	60	R3/38
10937/20	100	R3/38
10937/04	125 (4.oz)	R3/38
10937/05	175	R3/48
10937/06	250 (8 oz)	R3/48
10937/07	500	R3/58



Ointment jars

Without Caps

Clear glass, extra wide mouth

Product No	Capacity-ml	Cap Size
10936/08	15	R3/38
10936/00	30	R3/38
10936/01	60	R3/51
10936/02	120	R3/58
10936/03	250	R3/83
10936/04	500	R3/100



Amber Uncapped

Amber glass, extra wide mouth

Product No	Capacity-ml	Cap Size
10938/00	15	R3/28
10938/01	30	R3/33
10938/02	60	R3/38
10938/20	100	R3/38
10938/03	125 (4 oz)	R3/38
10938/07	175	R3/48
10938/04	250 (8 oz)	R3/48
10938/05	500	R3/58
10938/06	1000	R3/63



Capped

Clear glass, extra wide mouth and bakelite screw cap with PVDC faced liner

Product No	Capacity-ml	Cap Size
10078/08	15	R3/38
10078/00	30	R3/38
10078/01	60	R3/51
10078/02	120	R3/58
10078/03	250	R3/83
10078/04	500	R3/100



Clear Capped

Clear glass, wide mouth and bakelite screw cap, with PVDC faced liner

Product No	Capacity-ml	Cap Size
10076/01	15	R3/28
10076/02	30	R3/33
10076/03	60	R3/38
10076/04	125	R3/38
10076/05	175	R3/48
10076/06	250	R3/48
10076/07	500	R3/58



Kilner Jars

The well known Kilner preserving jar has a very wide mouth, and a resealable and replaceable metal lid. Ideal for solid samples, these bottles lend themselves to autoclaving due to the unique three part lid. Supplied singly in the capacities below, complete with lid.



Product No	Description	Capacity-ml
10953/01	Kilner Jar Complete	500
10953/02	Kilner Jar Complete	1000
10953/03	Sealing Ring, Pack of 12	
10953/04	Sealing Disc, Pack of 1	

Amber Capped

Amber glass, wide mouth and bakelite screw cap with PVDC faced liner

Product No	Capacity-ml	Cap Size
10077/00	15	R3/28
10077/01	30	R3/33
10077/02	60	R3/38
10077/03	125	R3/38
10077/07	175	R3/48
10077/04	250	R3/48
10077/05	500	R3/58
10077/06	1000	R3/63

