

Gas Jars

Simax clear glass gas jars with foot and ground flange.

Product No	Dia-mm	Height-mm
10503/01	50	150
10503/02	50	200
10503/03	50	250
10503/04	60	150
10503/05	60	300
10503/06	80	300



Discs for Gas Jars

Glass discs ground one side. Supplied singly

Product No	Dia-mm
10504/01	75
10504/02	115

Culture Bottles

Clear flint glass in the popular patterns indicated below.



With aluminium cap & rubber liner

Product No	Type	Capacity-ml	Cap Size	Pack Qty
10266/01	Bijou	7 (1/4 oz)	R3/20	288
10266/02	Bijou	14 (1/2 oz)	R3/20	288
10266/03	McCartney	28 (1 oz)	R3/22	144
10266/04	Universal	28 (1 oz)	R3/22	144

Without Caps

Product No	Type	Capacity-ml	Cap Size	Pack Qty
10083/01	Bijou	7 (1/4 oz)	R3/20	288
10083/02	Bijou	14 (1/2 oz)	R3/20	288
10083/03	McCartney	28 (1 oz)	R3/22	144
10083/04	Universal	28 (1 oz)	R3/22	144

Culture Bottles, Tubular

Thin walled, neutral borosilicate glass, tubular bottles, pre-capped with autoclavable polypropylene cap. Supplied in case quantities shown. Non-sterile. Discounts are available for bulk users.

Product No	Type	Capacity-ml	Pack Qty
10087/01	Universal	30	375
10087/02	Bijou	7	1715

Dreschel Bottles & Heads

Dreschel Bottles

Gas washing bottles for use with the Dreschel heads (supplied separately) below. All joints are 24/29.

Product No	Capacity-ml
10657/01	125
10657/02	250
10657/03	500



Plain Heads

Plain tube gas washing heads for use with above bottles.

Product No	Bottle Capacity-ml
10658/01	125
10658/02	250
10658/03	500



Adjustable Heads

Product No	Description
10659/01	Similar to 10658 above but with a screwthread joint on the receiving arm

Sintered Heads

Similar to Cat No 10658 above but available with a range of sinters at the end of the receiving arm to assist gas distribution. Suitable for use with bottle sizes indicated below.

Product No	Sinter Porosity	Bottle Capacity-ml
10660/01	0	125
10660/02	0	250
10660/03	0	500
10660/04	1	125
10660/05	1	250
10660/06	1	500
10660/07	2	250
10660/08	2	500
10660/09	3	250



Culture Bottle Racks

Autoclavable nylon coated wire bottle racks suitable for culture bottles containing agar slopes etc. Each rack has positions for two rows of six bottles of the sizes indicated below. Other rack sizes are available to order. Supplied singly.

Product No	Bottle Type
10437/01	Bijou
10437/02	Universal/McCartney