

Bottles and Jars

Dropping Bottles

Glass Dropping Bottles

Clear or Amber glass dropping bottles supplied complete with polystop stopper with integral teat and glass dropper.

Product No	Capacity-ml	Stopper.Size	Clear/Amber
10272/01	30	16/16	Clear
10272/02	50	16/16	Clear
10272/03	100	16/16	Clear
10272/21	50	16/16	Amber
10272/22	125	19/17	Amber

Density bottles

Soda glass, with drilled stoppers available adjusted (to BS733) or unadjusted.



Product No	Capacity-ml	Tolerance (+/- ml)	Adjusted Unadjusted
10567/01	10	0.5	A
10567/02	25	1.0	A
10567/03	50	1.0	A
10567/04	100	2.0	A
10567/05	10	0.5	U
10567/06	25	1.0	U
10567/07	50	1.0	U
10567/08	100	2.0	U

Medical Flat Bottles

Without Caps

Clear glass, screw neck. For caps, see accessory cap section

Product No	Capacity-ml	Cap Size
10081/01	100	R3/24
10081/02	150	R3/24
10081/03	200	R3/28
10081/04	300	R3/28
10081/05	500	R3/28



Capped

Clear glass, with rubber lined aluminium cap

Product No	Capacity-ml	Cap Size
10267/01	100	R3/24
10267/02	200	R3/28
10267/03	300	R3/28
10267/04	500	R3/28



Dispensing Bottles

Glass bottles fitted with Polystop stoppers, graduated droppers and teat.



Product No	Capacity/Dropper-ml	Clear/Amber	Stopper Size	Pack Qty
11034/01	100/1	Clear	16/16	10
11034/02	250/8	Clear	24/20	5
11034/03	500/10	Clear	24/20	5
11034/10	100/1	Amber	16/16	10
11034/11	250/8	Amber	24/20	5
11034/12	500/10	Amber	24/20	5

Weighing Bottles

Clear glass with ground in hollow stopper.

Tall Form

Product No	Height x Dia-mm	Pack Qty
10098/01	40 x 20	10
10098/02	50 x 25	10
10098/03	60 x 30	10
10098/04	80 x 40	5
10098/05	100 x 50	5



Squat Form

Product No	Height x Dia-mm	Pack Qty
10099/01	25 x 50	4
10099/02	30 x 60	4

Reagent Bottles, Blank Label

Laboratory bottles with plastic "polystop" stoppers, in clear or amber with a safety plastic coating and a blank label.

Product No	Capacity-ml	Clear/Amber
10376/01	250	Clear
10376/02	250	Amber
10376/03	500	Clear
10376/04	500	Amber