

Bottles and Jars

Reagent Bottles, Stoppered

Wide mouth, soda lime glass, with ground glass dustproof stopper.

Product No	Capacity-ml	Neck Size
10097/01	50	24/20
10097/02	100	29/22
10097/03	250	34/24
10097/04	500	45/27
10097/05	1000	60/46



Winchester Bottles

Narrow mouth, screw neck WITHOUT CAP (see later)

Clear Glass

Product No	Capacity	Cap Size
10968/01	30	R3/28
10968/02	60	R3/28
10968/03	90	R3/28
10968/04	100	R3/28
10968/05	125	R3/28
10968/06	150	R3/28
10968/07	180	R3/28
10968/08	200	R3/28
10968/09	250	R3/28
10968/10	300	R3/28
10968/11	500	R3/28
10968/12	1000	R3/28



Amber Glass

Product No	Capacity	Cap Size
10080/01	10	R3/20
10080/02	15	R3/20
10080/03	30	R3/20
10080/04	50	R3/24
10080/05	100	R3/24
10080/06	150	R3/24
10080/07	250	R3/28
10080/08	300	R3/28
10080/09	500	R6/31
10080/10	1000	R6/31
10080/12	2500	R6/31



Amber Wide Neck

Glass bottles with a 45mm neck (29mm I.D. for easy filling. Caps (10020 or 11050 range) must be purchased separately

Product No	Capacity
11049/01	500
11049/02	1000
11049/03	2500
11049/10*	2500

* PVC coated for added protection

Polypropylene caps for the 11049 range

Autoclavable to 140°C

Product No	Colour	Pack Qty
10020/03	Blue	10
10020/04	Yellow	10
10020/05	Green	10
10020/06	Grey	10
10020/07	White	10

11050/01* Blue 12

*Tamper-Evident Caps


Pouring Rings for the 11049 range

Autoclavable to 140°C. Supplied in packs of 10

Product No	Colour
10022/02	Blue
10022/03	Transparent

Specimen Tubes

Flat bottom soda glass tubes with polythene push plugs. Supplied in the dimensions and pack sizes indicated below.



Product No	Dia x Height-mm	Pack Qty
10279/01	38 x 10	250
10279/02	50 x 10	250
10279/03	38 x 12	100
10279/04	50 x 12	100
10279/05	75 x 12	100
10279/06	100 x 12	100
10279/07	50 x 16	100
10279/08	75 x 16	100
10279/09	50 x 18	100
10279/10	75 x 18	100
10279/11	50 x 25	100
10279/12	75 x 25	100
10279/13	100 x 25	100
10279/14	150 x 25	100