

Bottles and Jars

Bottles, Wide Mouth

Polyethylene, Round

Wide mouth bottle with DIN standardised screw caps, graduated. Supplied in pack quantities indicated below.



Product No	Capacity-ml	Pack Qty	Neck Size
10103/01	50	10	GL32
10103/02	100	10	GL32
10103/03	250	10	GL45
10103/04	500	10	GL45
10103/05	1000	5	GL63
10103/06	2000	1	GL63

Bottles, Extra Wide Mouth

HDPE

Moulded in high density natural polythene, squat form, with polypropylene screw cap.



Product No	Capacity-ml	Dia-mm	Height-mm	Neck ID-mm
10109/01	150	60	68	47
10109/02	285	78	92	57
10109/03	685	117	86	89
10109/04	1285	117	148	89

Security Bottles, Wide Mouth

TEFLON FEP, Round

High quality Nalgene leakproof bottles manufactured from Teflon FEP with a Tefzel ETFE screw cap. The Teflon construction means that these bottles can be used in the temperature range -100°C to +150°C, are autoclavable, and suitable for cryogenic work. Teflon is also highly chemically resistant and these bottles are therefore ideal for use with organic solvents.



Product No	Capacity-ml
10371/01	125
10371/02	250
10371/03	500
10371/04	1000

Polyethylene

A superior range of polyethylene bottles with a push-in insert under the cap to give a leakproof closure. The space between the cap and insert can be used for sample notes. Both bottle and cap have a small ring for security tagging and or sealing.



Product No	Capacity-ml	Pack Qty
10106/01	50	10
10106/02	100	10
10106/03	250	10
10106/04	500	10
10106/05	1000	5
10106/06	2000	1

Polypropylene, Round

High quality Nalgene leakproof bottles manufactured from autoclavable polypropylene. Supplied complete with screw cap and neck insert.



Product No	Capacity-ml	Pack Qty
10357/01	30	12
10357/02	60	12
10357/03	125	12
10357/04	250	12
10357/05	500	12
10357/06	1000	6