

Security Bottles, Narrow Mouth

Polyethylene

A superior range of polyethylene bottles, with a push plug insert under the cap which gives a leakproof closure. Both bottle and cap have a small hole for security sealing or tagging. Supplied complete with cap and insert.



Product No	Capacity-ml	Pack Qty
10105/01	50	10
10105/02	125	10
10105/03	250	10
10105/04	500	10
10105/05	000	5
10105/06	2000	1

Grey Polyethylene

Narrow neck bottles similar to Product No 10105 across except moulded in opaque grey polyethylene. Particularly suitable for transportation and storage of light sensitive samples.



Product No	Capacity-ml	Pack Qty
10107/01	50	10
10107/02	125	10
10107/03	250	10
10107/04	500	10
10107/05	1000	5
10107/06	2000	1

Bottles, Square

PVC

Leakproof, lightweight with polypropylene neck. Excellent for organic and inorganic acids, oxidising agents, aldehydes, amines, aliphatic hydrocarbons and mineral oils. Also useful for powders, granules and pastes.



Product No	Capacity-ml	Neck ID-mm	Pack Qty
11043/01	50	23	10
11043/02	100	31	10
11043/03	300	41	10
11043/04	500	54	10
11043/05	1000	69	5
11043/06	2000	69	1

Bottles, Square, Wide Mouth

Polypropylene

High quality Nalgene leakproof bottles manufactured from autoclavable polypropylene. The square shape of these bottles saves space, and is particularly suitable for long term storage or sample transportation. Supplied complete with cap and neck insert.



Product No	Capacity-ml	Pack Qty
10360/01	60	12
10360/02	175	12
10360/03	250	12
10360/04	500	12
10360/05	1000	6

