

Bottles and Jars

Bottles , Narrow Mouth

HDPE

High quality Nalgene leakproof bottles manufactured from almost unbreakable HDPE with a polypropylene screw cap. Supplied in the capacities and pack quantities indicated below. Note: 500ml size is available in two neck sizes.



Product No	Capacity-ml	Neck Size-mm	Pack Qty
10365/02	4	13	12
10365/03	8	20	12
10365/04	15	20	12
10365/05	30	20	12
10365/06	60	20	12
10365/07	125	24	12
10365/08	175	24	12
10365/09	250	24	12
10365/10	500	28	12
10365/11	500	38	12
10365/12	1000	38	6

HDPE, Amber

High quality Nalgene leakproof bottles manufactured from Amber HDPE with an Amber polypropylene screw cap. Complies with U.S. Pharmacopoeia XXII requirements for maximum light transmission. Supplied in the capacities and pack quantities indicated below.



Product No	Capacity-ml	Neck Size-mm	Pack Qty
10367/02	4	13	12
10367/03	8	20	12
10367/04	15	20	12
10367/05	30	20	12
10367/06	60	20	12
10367/07	125	24	12
10367/08	250	24	12
10367/09	500	28	12
10367/10	1000	38	6

LDPE

High quality Nalgene leakproof bottles manufactured from flexible L.D.P.E. with a polypropylene screw cap. Supplied in the capacities and pack quantities indicated below. Please note the 500ml size is available in two neck sizes.



Product No	Capacity-ml	Neck Size-mm	Pack Qty
10366/01	8	20	12
10366/02	15	20	12
10366/03	30	20	12
10366/04	60	20	12
10366/05	125	24	12
10366/06	250	24	12
10366/07	500	28	12
10366/08	500	38	12
10366/09	1000	38	6

Polyethylene

Graduated narrow mouth bottle with DIN standardised screw caps.



Product No	Capacity-ml	Pack Qty	Neck Size
10101/01	50	10	GL18
10101/02	100	10	GL18
10101/03	250	10	GL25
10101/04	500	10	GL25
10101/05	1000	5	GL32

Polypropylene

As Cat No 10101/xx series above, but polypropylene, autoclavable to 121°C.

Product No	Capacity-ml	Pack Qty	Neck Size
10102/01	50	10	GL18
10102/02	100	10	GL18
10102/03	250	10	GL25
10102/04	500	10	GL25
10102/05	1000	5	GL32

YORLAB