

## Bottles, Narrow Mouth

### Polypropylene

High quality Nalgene leakproof bottles manufactured from polypropylene, autoclavable to 121°C, and chemically resistant. Supplied complete with screw cap.



Product No	Capacity-ml	Neck Size-mm	Pack Qty
10368/02	4	13	12
10368/03	8	20	12
10368/04	15	20	12
10368/05	30	20	12
10368/06	60	20	12
10368/07	25	24	12
10368/08	250	24	12
10368/09	500	28	12
10368/10	1000	38	6

### Teflon FEP

High quality Nalgene leakproof bottles manufactured from Teflon FEP with a Tefzel ETFE screw cap. The Teflon construction means that these bottles can be used in the temperature range -100°C to +150°C, are autoclavable, and suitable for cryogenic work. Teflon is also highly chemically resistant and these bottles are therefore ideal for use with organic solvents.



Product No	Capacity-ml
10372/01	30
10372/02	60
10372/03	125
10372/04	250
10372/05	500
10372/06	1000
10372/07	2000

## Bottles, with Handle

### Polypropylene

White polypropylene bottles with integral handle and tamper evident cap. Particularly suitable for storage of chemical solutions.



Product No	Capacity-ml
10133/01	500
10133/02	1000
10133/03	2500
10133/04	5000

## Large Capacity Bottles

### HDPE

Narrow-mouth containers & caps, with carrying handle.

Product No	Capacity-Litres
11007/01	2.5
11007/02	5
11007/03	10



### HDPE

High quality Nalgene leakproof wide mouth bottles of a larger capacity manufactured from HDPE with a polypropylene screw cap. The neck I.D. is 88mm making filling easy.

Product No	Capacity-ml
10364/01	2000
10364/02	4000

## Jerricans

### HDPE

Rugged jerricans with a long spout for easy pouring, and a cap retaining strap.

Product No	Capacity-Litres
10432/01	5
10432/02	10
10432/03	20



## Autoclave Bottles

### Polypropylene

Wide-mouth, square 4-litre bottle with moulded hand-grips, and polypropylene cap. Excellent for storing and sterilising contaminated slides

Product No	11009/01
------------	----------

