

Aspirator Bottles, "Lowboys TM"

Nalgene "Lowboys" are manufactured from H.D.P.E., are leakproof, and because of their shape are suitable for storage and dispensing of liquids in a fridge. Can be stacked 3-high. Supplied complete with cap and tap.



Product No	Capacity (Litres)
10431/01	8
10431/02	15

Babcock Bottles

Glass babcock bottles with scale and writing patch. Supplied singly.

Product No	Description
10575/01	9g cream 60%
10575/02	18g cream 30%
10575/03	17.6ml milk 8%
10575/04	9g cream 50%



Bottles, B.O.D.

Bottles with a recessed stopper for biological oxygen demand determinations. Supplied singly in the capacities indicated below.

Product No	Capacity-ml
10400/01	250
10400/02	500



Bottles, Carboy

Larger volume robust H.D.P.E. bottles, narrow mouth with carrying handles. A push plug insert under the cap gives a superior seal. Supplied singly.

Narrow Mouth

Product No	Capacity (Litres)
10111/01	5
10111/02	10
10111/03	25
10111/04	50



Wide Mouth

Product No	Capacity (Litres)
10112/01	5
10112/02	10



Bottles, Carboy, Safety

Robust H.D.P.E. bottles with carrying handles and additional side neck for smooth easy pouring. Supplied singly.

Product No	Capacity (Litres)
10113/01	30
10113/02	50



Carboys, Polypropylene

Graduated polypropylene autoclavable carboys with handles and wide mouth leakproof screw closures. Supplied singly in the capacities indicated below.

100mm Closure

Product No	Capacity (Litres)
10426/01	10
10426/02	20



83B Closure

These carboys will accept filling/venting closures

Product No	Capacity (Litres)
10427/01	10
10427/02	20
10427/03	50



Carboys, Clear Polycarbonate

Nalgene leakproof polycarbonate "Clearboys™" are manufactured from tough but lightweight polycarbonate which is as clear as glass. With a leakproof 83B closure size these carboys will accept the filling/venting closures shown later. Autoclavable and graduated. Supplied singly in the capacities indicated below.

Product No	Capacity (litres)
10428/01	10
10428/02	20



Disposable Bijou, Sterile

Sterile polystyrene 7ml bijou containers, with screw caps providing a useful write-on area. Supplied with or without a plain white label on the bottle in packs of 700.

Product No	Label
11033/01	No
11033/02	Yes

