

Caps for Low Cost Bottles

For 10130 series bottles, supplied in packs of 1000

Product No	Tamper evident	Colour
10132/05	Yes	Blue

Filling/Venting Closures

All-plastic closures for aseptic liquid transfer of media, biological reagents, pure water and chemicals to and from Nalgene bottles and carboys with a 53mm (53B) or 83mm (83B) necks. Autoclavable.



Product.No	Neck Size	Tubing Size-ID
10433/01	53B	6mm
10433/02	83B	12mm
10433/03	83B	6mm

Jars

HDPE

Cylindrical jars with caps & inserts. The large neck makes them particularly suitable for powders, solids and semi-solids.



Product No	Capacity-ml	Neck ID-mm	Pack Qty
11055/01	70	36	10
11055/02	120	36	10
11055/03	250	50	10
11055/04	500	70	10
11055/05	1000	85	5
11055/06	1500	87	1
11055/07	2000	87	1

Jars, Wide Mouth

Polypropylene

Translucent polypropylene with excellent chemical resistance and repeatedly autoclavable. With polypropylene cap.



Product No	Capacity-ml	Closure size-mm	Pack Qty
11042/01	15	38	12
11042/02	30	43	12
11042/03	60	53	12
11042/04	125	70	12
11042/05	250	70	6
11042/06	500	120	6
11042/07	1000	120	6

Seals With Aperture

Silicone rubber sealing rings, P.T.F.E. faced allow GL thread aperture caps to have a variety of implements permanently sealed in place at the top of a bottle or any item carrying a GL thread.

Examples are Kipps type dispensers, thermometers, stopcocks etc. The pressure of screwing down the cap squeezes the silicone to grip the item. Supplied in packs of 10 for the cap sizes and implement diameters shown below.



Product No	Cap Size	Implement.Dia
10516/01	GL14	5.5 to 6.5
10516/02	GL18	5.5 to 6.5
10516/03	GL18	7.5 to 9.0
10516/04	GL18	9.0 to 11.0
10516/05	GL25	7.5 to 9.0
10516/06	GL25	9.0 to 11.0
10516/07	GL25	11.0 to 13.0
10516/08	GL32	9.0 to 11.0
10516/09	GL32	11.0 to 13.0
10516/10	GL32	13.0 to 15.0
10516/11	GL32	15.0 to 17.0
10516/12	GL32	17.0 to 19.0
10516/13	GL45	25.0 to 27.0
10516/14	GL45	31.0 to 33.0

Seals, Pierceable

Pierceable silicone rubber seals for use with aperture caps. Supplied in packs of 10.

Product No	Cap Size
10091/01	GL14
10091/02	GL18
10091/03	GL25
10091/04	GL32
10091/05	GL45

Seals, PTFE Faced

P.T.F.E. faced silicone rubber seals for use with aperture caps. Will withstand temperatures of 180°C. Supplied in packs of 10.

Product No	Cap Size
10092/01	GL14
10092/02	GL18
10092/03	GL25
10092/04	GL32
10092/05	GL45