


Burettes, Auto Zero, Plastic


Acrylic bodied burettes. Supplied complete with an L.D.P.E. squeeze bottle and the necessary tubing for simple burette filling.

Product No	Capacity-ml
10930/01	10
10930/02	25
10930/03	50
10930/04	100



Burette Clamps

A polypropylene assembly that will fit support rods with diameters between 8-14mm, and holding one or two burettes, as indicated. Spring loaded grips hold the burette securely while allowing easy removal. Graduations and meniscus remain visible. Suitable retort rods and bases may be found in the Stands & Supports section, if required. Not supplied with support rod.



Product No	No of Positions
10835/01	1
10835/02	2

Bunsen Tubes

Neoprene tubes with ends reinforced for repeated connection/disconnection of burners. 7-8mm bore.

Product No	Length-mm
10340/01	450
10340/02	600
10340/03	750
10340/04	900



Gas Lighter

Piezo-electric operation, requiring no battery

Product No 11006/01



Spirit Burners

Glass spirit burners with fitted wick in ceramic burner. Supplied complete with snuffer cap.

Product No	Capacity-ml
10456/01	50
10456/02	100
10456/03	spare wick 10m



Gauzes plain

Suitable for use with tripods, plain 20 mesh iron wire gauze. Supplied singly in the sizes indicated below.

Product No	Size-mm
10346/01	100 x 100
10346/02	125 x 125
10346/03	150 x 150



Gauzes, Ceramic Centre

Iron wire gauzes with ceramic centres supplied singly in the sizes indicated below.

Product No	Size-mm
10347/01	100 x 100
10347/02	125 x 125
10347/03	150 x 150
10347/04	200 x 200



Pipe clay, Triangles

Tubes of pipe clay on iron wire. Supplied singly in the side lengths indicated below.

Product No	Length-mm
10348/01	40
10348/02	50
10348/03	65
10348/04	75
10348/05	90

