

Ranger™ High resolution - Fast, high precision balance in a rugged industrial case

Applications

- Weighing (4 units)
- Mains operation or optional battery operation

Features

- RS232 interface
- Backlit LCD display
- Net/gross/tare function
- Auto-timer shut-off
- Levelling feet and bubble
- Rugged cast aluminium and steel construction
- Refer to later pages for the correct calibration weights

Options

- Internal rechargeable battery
- EC Type Approval - second stage
- Second RS232

Product no	Capacity	Factory readability	Approved Readability	Platform
12306/01	1.5kg	0.01g	0.1g	165 x 165
12306/02	3kg	0.01g	0.1g	165 x 165
12306/03	6kg	0.1g	1g	200 x 240
12306/04	15kg	0.1g	1g	240 x 350
12306/05	35kg	0.1g	1g	240 x 350



Options

Must be ordered at the same time as the scale

Product no	Description	
12206/01	2nd stage verification*	Add "1" to the end of the Model Number (eg RC3RD1)
12206/02	Additional RS232	Add "2" to the end of the Model Number (eg RC3RD2)
12206/03	Internal Battery	Add "3" to the end of the Model Number (eg RC3RD3)

* Postcode of final destination must be supplied at time of order

Accessories

Order no	Description
12206/10	Anti theft device
12206/11	RS232 cable for 2nd scale 1.8m
12206/12	Carry case for small platform models
12206/13	Carry case for large platform models
12206/14	Ranger working cover for small platform
12206/15	Ranger working cover for large platform
12206/16	Replacement AC mains adapter
12206/17	Printer SF42 (includes interface cable)
12206/18	PC Interface cable - 9 pin

