

Burners and Butryometers

Tripods Triangular

Nickel plated mild steel tripods with PVC feet.

Product No	Triangle side-mm	Height-mm
10836/01	125	200
10836/02	150	200



Bunsen Burners

Overall height approximately 124 mm. dimension given is the outside diameter of the burner tube.

With air regulator

	o.d. mm	for gas type
12290/01	11	Natural
12290/02	13	Natural
12290/03	11	LPG
12290/04	13	LPG

With control stopcock and pilot flame

	o.d. mm	for gas type
12290/05	13	Natural
12290/06	13	LPG



Meker Burners, Meker pattern

With metal grid to ensure an evenly distributed flame. Overall height approximately 140mm. Dimension is the grid diameter.

	o.d. mm	for gas type
12290/10	30	Natural
12290/11	43	Natural
12290/12	30	LPG
12290/13	43	LPG



Butryometers, B.S.

Manufactured to BS696 suitable for testing dairy products by the Gerber method. For 10.94ml pipettes.

Product No	Description
10293/01	Cheese 40% 3g
10293/02	Cream 70% 5g
10293/03	Skimmed Milk 4% 10.94ml
10293/04	Milk 8% 10.94ml
10293/05	Skimmed Milk Double Qty 0.5%

Butryometers, I.S.

Manufactured to Irish Standard specification (IS66 & 68) suitable for testing dairy products by the Gerber method. For 10.8ml pipettes.

Product No	Description
10294/01	Cheese 40% 3g
10294/02	Cream 70% 5g
10294/03	Milk 7% 10.8ml
10294/04	Milk 8% 10.8ml

Butryometer Lockstoppers

Product No	Description
11037/01	With brass rim, solid glass locking ball and long white rubbers

Butryometer Brushes

Product No	Description
10379/01	Brushes specially shaped for cleaning butryometers. Supplied singly

