

## Microcentrifuge tubes and more z 233 M-2 / Z 233 MK-2



### Features

- Brushless drive
- Microprocessor control
- Imbalance switch off
- Lid lock

### Digital Display for:

- Speed
- Force
- Running time
- Temperature
- Key for short runs
- Acceleration/deceleration profile
- Continuous run function
- Standstill cooling
- Preselection by potentiometers
- Diagnosis programme

### Specifications Z 233 M-2

Max. volume:	44 x 1.5/2.0ml
Max speed:	15000 rpm
Max. RCF:	21380 xg

### Dimensions

w x h x d (cm):	28 x 26 x 37
Weight:	17.0 kg

**Product no** 12370/30

### Specifications Z 233 MK-2

#### Refrigerated

Max. volume:	6 x 30ml
Max speed:	15000 rpm
Max. RCF:	21380 xg

### Dimensions

w x h x d (cm):	32 x 29 x 54
Weight:	39.0 kg

**Product no** 12370/31

### Accessories

Product no	Description			
12370/32	Angle rotor	18 x 1.5/2.0ml	14000 rpm	15996 xg
12370/33	Angle rotor	24 x 1.5/2.0ml	15000 rpm	21380 xg
12370/34	Angle rotor	44 x 1.5/2.0ml	14000 rpm	18400 xg
12370/35	Angle rotor	64 x 0.4ml	14000 rpm	17960 xg
12370/36	Angle rotor	4 x 1.5/0.2ml strips		15339 xg
12370/37	Hematocrit rotor	24 x capill. 1.3 x 75mm	14811 xg	
12370/38	Angle rotor	6 x 30ml	12000rpm	15938 xg

can only be used in the refrigerated Z 233 MK-2