

Centrifuges

Universal 4 x 100ml



Features

- Brushless drive
- Microprocessor control
- Imbalance switch off
- Lid lock
- Automatic rotor recognition

Digital Display for:

- Speed
- Force
- Running time
- Temperature
- Key for short runs
- Acceleration/deceleration profile
- Continuous run function
- Standstill cooling
- Preselection by potentiometers
- Diagnosis programme

Specifications Z 300

Max. volume:	4 x 100ml
Max speed:	13500 rpm
Max. RCF:	17320 xg
Dimensions	
w x h x d (cm):	35 x 33 x 37
Weight:	17.0 kg
Product no	12370/40

Specifications Z 300 K

Refrigerated

Max. volume:	4 x 100ml
Max speed:	13500 rpm
Max. RCF:	17320 xg
Dimensions	
w x h x d (cm):	42 x 33 x 66
Weight:	55 kg
Product no	12370/41

Accessories

Product no	Description			
12370/42	Swing out rotor	4 x 100ml	4500 rpm	3350xg
12370/43		8 x 15ml Falcon		
12370/44		36 x 15ml		
12370/45	Swing out rotor	2 x 3 Microtitre plates	4500 rpm	2762 xg
12370/46	Angle rotor	6 x 50ml	6000 rpm	3904 xg
12370/47	Angle rotor	12 x 15ml	6000 rpm	4185 xg
12370/48	Angle rotor	24 x 15ml	3500 rpm	2000 xg
12370/49	Angle rotor	44 x 1.5/.2.0ml	13500 rpm	17320 xg
12370/42	Hematocrit rotor	24 x capill. 1.3 x 75mm		14811 xg

YOR

email sales@yorlab.co.uk